

YOUR MEDICAL HARDWARE STORE  
708-369-3675 - **24 HOUR SERVICE**

888-998-VMED (8633)

**\* RENTALS \* REPAIRS \***  
**\* SALES \* SHUTDOWNS \***



# Value Med



SERVICES



RESPIRATORY



SUCTION



TECH SUPPORT



## PLEASE CONTACT US AT:

EMAIL: ValueMed.sales@gmail.com  
WEBSITE: www.valuemed.co  
TOLL FREE: 888-998-VMED (8633)  
FAX: 888-676-0711

NOTES:

---



---



---



---



---



**Kevin McDermott**  
 (708) 369-3675  
 EMAIL: KMCDMED@aol.com  
 FAX: 888-676-0711

**Value Med**

**INDEX**

Air Oxygen Blenders.....	12-13	Outlet Accessories.....	52, 53
Adapters .....	6-7	Oxygen Sensors / Fuel Cells.....	16
Ball Valves / Check Valves.....	47	Oxygen Analyzers / Monitors.....	17
Bulb/Lamp, Replacements.....	20	Percussors.....	15
Cufflators.....	15	Pressure & Vacuum Gauges.....	52-53
Cylinder Carts / Stands.....	23	Rail Systems.....	39
Diss Body Adapters.....	44-45	Respiratory Accessories.....	21
Flowmeters:		Rolling Stands.....	22
Standard.....	8-9		
Specialty.....	32-33		
High Pressure Hoses.....	40-41		
Insufflator / Pigtails.....	38		
Laryngoscopes.....	18-19		
NIF Gauges.....	14		
Nitrogen Accessories.....	36-37		
NPT Fittings.....	46		



Regulators:

- Carbon Dioxide.....35
- Oxygen.....34
- Medical Air.....34
- Nitrous Oxide.....34
- Heli-Ox.....34
- Nitrogen.....36

Spirometers.....15

Suction Accessories.....31

Suction Canisters / Rings.....28-29

Suction Pumps.....30

Suction Regulators:

- Boehringer.....24
- Chemetron.....25
- Precision Medical.....25
- Ohio / Ohmeda.....27
- Amvex.....26

Wall Outlets:

- Puritan Bennett /  
Beacon Medaes.....49
- Chemetron.....51
- Medstar.....51
- Ohio.....50

Wye Assemblies.....7

**ABOUT VALUEMED**

Welcome to Valuedmed,  
 Valuedmed business encompasses hospital sales, biomedical services, onsite medical gas verification, hospital rental equipment, as well as being the largest medical gas and vacuum shutdown supplier in the midwest.  
 Valuedmed’s sales and technical specialists are available to respond to the unique situations involving the ever changing needs of hospitals.

***Valuedmed Professional Services***

- \*Custom High Pressure Hose Assemblies.....40,41
- \*Onsite Medical Gas Analysis.....48
- \*Outlet Verification.....48
- \*Equipment Repair.....4
- \*Equipment Rental.....5
- \*Shutdown Services.....55



Thank You For Allowing Us to Be Your *VALUED* Partner In Healthcare

**EQUIPMENT REPAIR SERVICES**

**\*\*CALL US TODAY FOR EQUIPMENT REPAIR PICK-UP\*\***

**Medical Gas  
Regulators & Flowmeters**



Oxygen  
Medical Air  
Heli-Ox



Carbon Dioxide  
Nitrogen  
Nitrous Oxide

**Suction Regulators**



**Repair or Upgrade Your  
Overflow Safety Trap**

**\$39.50**

**Call For  
Estimates**



\* Chemetron \* Ohmeda \* \* Precision Medical \* Boehringer \*

**Air/Oxygen Blenders**

Overhauls  
Water Traps  
Hoses  
Brackets



Bird - Bio Med Devices - Precision Medical  
*Recommended Overhall every 2 years*

**FULL LINE OF CUSTOM**

**HOSES AVAILABLE**

\*High Pressure Hoses \*Hose  
Assemblies\*Surgery\*

**\*Anesthesia\*Respiratory\***

**Thank You Kevin McDermott**

**(708)369-3675**

**KMCDMED@AOL.COM**



**Oxygen Analyzers & Monitors  
Replacement Sensors**



MiniOx \* Hudson \*Maxtec

**Ventilator Support Arms**



**Horizon Gas Block Sales and Service**



**Laryngoscope  
Blades & Handles**

\* Replacement\*  
Lamps



**Cufflators**



Manometers

**Respirometers**



Wright -Boehringer-Haloscale

**Percussors**



General Physiotherapy  
Meds Systems





**YOUR MEDICAL HARDWARE STORE**  
**708-369-3675 - 24 HOUR SERVICE**  
**888-998-VMED (8633)**

## RESPIRATORY EQUIPMENT

**\* RENTALS \* REPAIRS \* SERVICE \* SHUTDOWNS \***



**VENTILATORS**  
 \*SERVO I \*840



**NON-INVASIVE VENTILATORS**  
 V60 BIPAP



BIPAP VISION

## SHUTDOWN SERVICES

### VAPOTHERM

PRECISION  
 FLOW WITH  
 IV STAND



### PATIENT CIRCUITS

PF-E-DPC LOW FLOW

PF-E-DPC HIGH FLOW

*LET US HELP WITH YOUR EMERGENCY AND PLANNED SHUTDOWNS*

### COMPRESSORS AND SUCTION PUMPS



### CYLINDER CARTS



*DUAL & SINGLE  
 CARTS AVAILABLE*

### REGULATORS



### COUGH ASSIST



### SHUTDOWN EQUIPMENT

- \* REGULATORS/ PIGTAILS
- \* CYLINDER CARTS
- \* HOSES
- \* AIR COMPRESSORS
- \* SUCTION PUMPS
- STAND & TABLE TOP

**CHEMETRON QUICK CONNECT ADAPTERS**



PART #	GAS	PRICE
401-55	O2	\$29.00
411-55	AIR	\$29.50
421-55	N2O	\$46.00
431-55	VAC	\$36.00
451-55	EVAC	\$89.00
481-55	CO2	

PART #	GAS	PRICE
401-18	O2	\$24.50
411-18	AIR	\$28.50
421-18	N2O	\$39.50
431-18	VAC	\$32.50
451-18	EVAC	\$89.00
481-18	CO2	



**PURITAN BENNETT QUICK CONNECT ADAPTERS**




PART #	GAS	PRICE
201-55	O2	\$26.00
211-55	AIR	\$29.00
221-55	N2O	\$46.00
231-55	VAC	\$36.00
251-55	EVAC	\$79.00
281-55	CO2	

PART #	GAS	PRICE
201-18	O2	\$15.50
211-18	AIR	\$18.50
221-18	N2O	\$39.50
231-18	VAC	\$32.50
251-18	EVAC	\$59.00
281-18	CO2	




**OHIO QUICK CONNECT ADAPTERS**



PART #	GAS	PRICE
501-55	O2	\$28.50
511-55	AIR	\$29.50
521-55	N2O	\$37.50
531-55	VAC	\$36.50
551-55	EVAC	\$69.00
581-55	CO2	\$89.00

PART #	GAS	PRICE
501-18	O2	\$27.00
511-18	AIR	\$27.50
521-18	N2O	\$39.00
531-18	VAC	\$31.50
551-18	EVAC	\$69.00
581-18	CO2	\$69.00



**MEDSTAR QUICK CONNECT ADAPTERS**



PART #	GAS	PRICE
302-55	O2	\$64.00
312-55	AIR	\$64.00
322-55	N2O	\$68.00
332-55	VAC	\$68.00
354-55	EVAC	\$89.00
355-55	CO2	\$89.00

PART #	GAS	PRICE
302-18	O2	\$49.50
312-18	AIR	\$49.50
322-18	N2O	\$64.00
332-18	VAC	\$55.00
354-18	EVAC	\$89.00
355-18	CO2	\$89.00



**OXEQUIP QUICK CONNECT ADAPTERS**



PART #	GAS	PRICE
301-55	O2	\$49.00
311-55	AIR	\$49.00
321-55	N2O	\$59.00
331-55	VAC	\$54.00

PART #	GAS	PRICE
301-18	O2	\$24.00
311-18	AIR	\$26.00
321-18	N2O	\$38.00
331-18	VAC	\$29.00





**Value Med**

**TO PLACE AN ORDER**

**FAX: 888-676-0711**

**EMAIL: [valuemed.sales@gmail.com](mailto:valuemed.sales@gmail.com)**

**WYE ASSEMBLIES**

**OHIO WYE ASSEMBLY**



**OXYGEN**

22-501-500		\$144.00
------------	--	----------

**AIR**

22-511-510		\$149.00
------------	--	----------

**VACUUM**

22-531-530		\$169.00
------------	--	----------

**CHEMETRON WYE ASSEMBLY**



**OXYGEN**

22-401-400		\$144.00
------------	--	----------

**AIR**

22-411-410		\$149.00
------------	--	----------

**VACUUM**

22-431-430		\$169.00
------------	--	----------

**PURITAN BENNETT WYE ASSEMBLY**



**OXYGEN**

22-201-200		\$144.00
------------	--	----------

**AIR**

22-211-210		\$149.00
------------	--	----------

**VACUUM**

22-231-230		\$169.00
------------	--	----------

**MEDSTAR WYE ASSEMBLY**



**OXYGEN**

22-302-300		\$199.00
------------	--	----------

**AIR**

22-312-310		\$199.00
------------	--	----------

**VACUUM**

22-332-330		\$199.00
------------	--	----------

**DISS WYE ASSEMBLY**



22-102-100	OXYGEN	\$144.00
------------	--------	----------

22-112-110	AIR	\$149.00
------------	-----	----------

22-132-130	VACUUM	\$169.00
------------	--------	----------

132EMG	Green	100/Pkg
0137	Yellow	50/Pkg
0135	White	50/Pkg
0135BLU	Blue	50/Pkg
0132BKW	Black/White	50/Pkg
0132BLK	Black	50/Pkg
0138	Clear	50/Pkg

**SWIVEL STYLE TUBING NIPPLE**



**TO ORDER: Phone: 888-998-8633 / Fax: 888-676-0711 / Email: [Valuemed.Sales@Gmail.com](mailto:Valuemed.Sales@Gmail.com)**















**Value Med**

**YOUR MEDICAL HARDWARE STORE**

**24 HOUR SERVICE**

**888-998-VMED (8633)**

## STANDARD AIR / OXYGEN FLOWMETERS

MODEL						
	0-15 LPM STANDARD O2 FLOWMETER	0-15 LPM W/ 50 PSI POWER TAKE OFF	0-15 LPM STANDARD AIR FLOWMETER	0-3 LPM PED/NEO OXYGEN *	0-1 LPM PED/NEO OXYGEN**	200 CC PED/NEO OXYGEN***
 OHIO	15-501	OXYGEN 15-501PT AIR 15-510PT	15-510	3-501	1-501	200-501
 PURITAN BENNETT	15-201	OXYGEN 15-201PT AIR 15-211PT	15-211	3-201	1-201	200-201
 CHEMETRON	15-401	OXYGEN 15-401PT AIR 15-411PT	15-411	3-401	1-401	200-401
 OXEQUIP	15-301	OXYGEN 15-301PT AIR 15-311PT	15-311	3-301	1-301	200-301
 MEDSTAR	15-302	OXYGEN 15-302PT AIR 15-312PT	15-312	3-302	1-302	200-302
 DISS	15-102	OXYGEN 15-102PT AIR 15-112PT	15-112	3-102	1-102	200-102
1/8" NPT FEMALE	102-18F	OXYGEN 102-18FPT AIR 15-112-18FPT	15-112-18	3-102-18	1-102-18	200-102-18

\* 0-3 LPM 1/8 LPM UP TO 1 LPM 1/4 LPM FROM 1-3 \*\* 0-1 LPM 1/10 LPM INCREMENTS \*\*\* 0-200 cc -25cc INCREMENTS

















**Value Med**

**TO PLACE AN ORDER**

**FAX: 888-676-0711**

**EMAIL: [valuemed.sales@gmail.com](mailto:valuemed.sales@gmail.com)**

## STANDARD AIR / OXYGEN FLOWMETERS


<b>MODEL</b>	<b>OXYGEN/AIR FLOWMETER WYE</b>	<b>OXYGEN/AIR FLOWMETER WYE W/ PTO 50 PSI</b>		<b>SECURITY CABLE IDEAL FOR FLOWMETERS</b>
	 OXYGEN / AIR FLOWMETER WYE	 OXYGEN / AIR FLOWMETER WYE WITH POWER TAKEOFF	<b>VENTILATOR / OUTLET SPACE SAVER</b>	 15-167   \$3.75 Each
 OHIO	OXYGEN 22-15-501 AIR 22-15-510	OXYGEN 22-15-501PT AIR 22-15-510PT	T15-501	<b>FLOW SELECTOR</b>  PM1000
 PURITAN BENNETT	OXYGEN 22-15-201 AIR 22-15-211	OXYGEN 22-15-201PT AIR 22-15-211PT	T15-201	
 CHEMETRON	OXYGEN 22-15-401 AIR 22-15-411	OXYGEN 22-15-401PT AIR 22-15-411PT	T15-401	
 OXEQUIP	OXYGEN 22-15-301 AIR 22-15-311	OXYGEN 22-15-301PT AIR 22-15-311PT	T15-301	
 MEDSTAR	OXYGEN 22-15-302 AIR 22-15-312	OXYGEN 22-15-302PT AIR 22-15-312PT	T15-302	
 DISS	OXYGEN 22-15-102 AIR 22-15-112	OXYGEN 22-15-102PT AIR 22-15-112PT	T15-100	
1/8" NPT	OXYGEN 22-102-18F AIR 22-112-18F	OXYGEN 22-102-18FPT AIR 22-112-18FPT	N/A	<b>POWER TAKE OFF 50 PSI</b>  OXYGEN P-24 AIR P-16


**\*\*TRADE IN BROKEN FLOWMETERS FOR ONLY \$22.50\*\***

**TO ORDER: Phone: 888-998-8633 / Fax: 888-676-0711 / Email: [Valuemed.Sales@Gmail.com](mailto:Valuemed.Sales@Gmail.com)**

**“E” & “H” TANK**


**50 PSI PRESET REGULATORS PURITAN BENNETT / CHEMETRON**


50 PSI PRESET MODEL WITH PB FEMALE			
			
50-870-200	“E” Tank, O2	870 CGA	Yoke
50-540-200	“H” Tank, O2	540 CGA	Nut
50-950-210	“E” Tank, Air	950 CGA	Yoke
50-346-210	“H” Tank, Air	346 CGA	Nut
50-910-220	“E” Tank, N2O	950 CGA	Yoke
50-236-220	“H” Tank, N2O	346 CGA	Nut

50 PSI PRESET MODEL WITH CHEMETRON FEMALE			
			
50-870-400	“E” Tank, O2	870 CGA	Yoke
50-540-400	“H” Tank, O2	540 CGA	Nut
50-950-410	“E” Tank, Air	950 CGA	Yoke
50-346-410	“H” Tank, Air	346 CGA	Nut
50-910-420	“E” Tank, N2O	950 CGA	Yoke
50-236-420	“H” Tank, N2O	346 CGA	Nut

**“E” & “H” TANK**

**50 PSI PRESET REGULATORS OHIO / MEDSTAR**

50 PSI PRESET MODEL WITH OHIO FEMALE			
			
50-870-500	“E” Tank, O2	870 CGA	Yoke
50-540-500	“H” Tank, O2	540 CGA	Nut
50-950-510	“E” Tank, Air	950 CGA	Yoke
50-346-510	“H” Tank, Air	346 CGA	Nut
50-910-520	“E” Tank, N2O	950 CGA	Yoke
50-236-520	“H” Tank, N2O	346 CGA	Nut

50 PSI PRESET MODEL WITH MEDSTAR FEMALE			
			
50-870-300	“E” Tank, O2	870 CGA	Yoke
50-540-300	“H” Tank, O2	540 CGA	Nut
50-950-310	“E” Tank, Air	950 CGA	Yoke
50-346-310	“H” Tank, Air	346 CGA	Nut
50-910-320	“E” Tank, N2O	950 CGA	Yoke
50-236-320	“H” Tank, N2O	346 CGA	Nut




**Value Med**


**TO PLACE AN ORDER**


**FAX: 888-676-0711**



**EMAIL: [valuemed.sales@gmail.com](mailto:valuemed.sales@gmail.com)**


**OXYGEN / AIR / NITROUS OXIDE REGULATORS**

50 PSI PRESET REGULATORS		
		
<b>OXYGEN</b>		<b>\$159.00</b>
<b>50-870</b>	Small Cylinder E Yoke	870 CGA
<b>50-540</b>	Large Cylinder Nut	540 CGA
<b>AIR</b>		<b>\$179.00</b>
<b>50-950-08</b>	Small Cylinder E Yoke	950 CGA
<b>50-346-08</b>	Large Cylinder Nut	346 CGA
<b>NITROUS OXIDE</b>		<b>\$289.00</b>
<b>50-910</b>	Small Cylinder E Yoke	910 CGA
<b>50-326</b>	Large Cylinder Nut	326 CGA

0-100 PSI ADJUSTABLE REGULATORS		
		
<b>OXYGEN</b>		<b>\$199.00</b>
<b>100-870</b>	Small Cylinder E Yoke	870 CGA
<b>100-540</b>	Large Cylinder Nut	540 CGA
<b>AIR</b>		<b>\$249.00</b>
<b>100-950</b>	Small Cylinder E Yoke	950 CGA
<b>100-346</b>	Large Cylinder Nut	346 CGA
<b>NITROUS OXIDE</b>		<b>\$289.00</b>
<b>100-910</b>	Small Cylinder E Yoke	910 CGA
<b>100-326</b>	Large Cylinder Nut	326 CGA

OXYGEN E-REGULATOR	
	
<b>PART#</b>	<b>GAS</b>
15-870	Standard 15 LPM Oxygen
25-870PT	25 Oxygen W/ 50psi PTO

Oxygen Regulator with Gauge	OXYGEN MRI REGULATOR
	
Oxygen Regulator	MRI Regulator / Flowmeter

"E" CYLINDER WRENCH & SECURITY CABLE		
		
<b>BMCW-3P</b>	Wrench Only	\$4.90
<b>BMCW-2BC</b>	Wrench & Cable	\$7.85

PART #	DESCRIPTION	PRICE
15-870-DI	"E" Cylinder, gauge, 15 LPM	\$129.00
15-540-DI	"H" Cylinder, gauge, 15 LPM	\$129.00
15-870-15	"E" Cylinder/15 LPM Flowmeter/Yoke	\$195.00
15-540-15	"H" Cylinder/15 LPM Flowmeter/Nut	\$189.00
*MRI Regulators include 0-15 LPM Flowmeters*		
15MRI-870-DI	Non-Ferrous "E" Regulator, 870 Yoke	\$195.00
15MRI-540-DI	Non-Ferrous "H" Regulator, 540 Nut	\$195.00

## AIR OXYGEN BLENDERS



**BIOMED DEVICES**



**TENACORE**



**PRECISION  
MEDICAL**



**BLENDER BUDDY**

Many Options Available  
Call for Details & Pricing

### \*\*BLENDER OVERHAULS\*\*

#### OVERHAULS INCLUDE

- \*Ultrasonic Cleaning
- \*Replacement Of Necessary Parts
- \*Testing \*\*Full PM Kit
- \*Calibration To Manufacturers Specs
- \*Bird
- \*BioMed Devices
- \*Tenacore
- \*Precision
- \*Sechrist

**Overhauls are Recommended Every 2 Years**

#### STANDARD AIR / OXYGEN BLENDER

#BM-700-01	\$799.00
------------	----------

#### Pediatric Low Flow Air / Oxygen Blender

#BM-700-03	\$839.00
------------	----------

#### MRI Air / Oxygen Blenders

#BM-700-1MRI	\$1,079.00
--------------	------------


#### REFURBISHED BLENDERS AVAILABLE

#### Blender Accessories

- BM-108-007 10 FT Oxygen Hose, DISS Male x DISS Female
- BM-208-007 10 FT Air Hose, DISS Female x DISS Female
- BM-300-007 Elbowed 0-15 LPM Oxygen Flowmeter
- BM-400-007 Elbowed Pediatric 0-3 LPM Oxygen Flowmeter
- BM-500-007 Elbowed 0-1 LPM Oxygen Flowmeter
- BM-570-007 Elbowed 0-70 LPM Oxygen Flowmeter
- BM-600-007 Elbowed 15 LPM Oxygen Flowmeter Wye
- AFA11629 Air Inlet Filter/Water Trap
- 1059-007 Stand with Wheels

#### Brackets

- BR-800-007 Pole Mount Bracket
- BR-825-007 Female Dovetail Bracket
- BR-850-007 Wall Mount Bracket
- BR-901-007 Hill Rom / Fairfield Blender Bracket
- BR-903-007 Suction Slide Bracket to Blender
- BR-904-007 Vertical Rail Bracket, GE Incubator



**GE  
VERTICAL  
RAIL BRACKET**

BR-904-007



**ELBOWED  
FLOWMETER**

BM-300-007



**STAND WITH  
WHEELS**

1059

#### POLE MOUNT BRACKET




BR-800-007

#### WATER TRAP



AFA116-2-9



**VAC BRACKET  
TO BLENDER**

BR-903-007



**HILL ROM  
BLENDER  
BRACKET**

BR-901-007



**WALL MOUNT  
BRACKET**

BR-850-007





**MAXFLO-2**  
**MF-HIGH FLO \$599**  
**MF-LO FLO \$599**



**MAX-BLEND W/DISS**  
**MB-HIGH FLO.....\$1599.00**  
**MB-LO FLO.....\$1799.00**



**MAX-VENTURI WITH SENSOR**  
**MV-S.....\$999.00**



**SECHRIST**  
**3500 CP-G**

**OXYGEN CONCENTRATORS**



**INVACARE**  
 10L Platinum 10  
 W/Sensor  
**IRC-10LX-O2**  
 5L PerfectO2 W/Sensor  
**IRC-5PW**




**DAMPENS NOISE AND MINIMIZES PATIENT STIMULATION**

**HIGH FLOW SILENCER**  
**FITS 15, 22MM HOSE & FLOWMETERS**  
**\$97.00 each.....219P50**

**Precision Endotracheal Tube Anti-Slippage Clamp**

**\* Modified Precision Cord Clamp \***



**\*CLAMP SOLD SEPARATE FOR MODIFIED CLAMP**

Part ET-20	ET Tube Size 2.0 MM	\$3.99 Each
Part ET-25	ET Tube Size 2.5 MM	\$3.99 Each
Part ET-30	ET Tube Size 3.0 MM	\$3.99 Each
Part ET-35	ET Tube Size 3.5 MM	\$3.99 Each
Part ET-40	ET Tube Size 4.0 MM	\$3.99 Each
Part ET-45	ET Tube Size 4.5 MM	\$3.99 Each
Part M100	Clamp Only	\$1.55 Each

**\*CLAMP SOLD SEPARATE FOR MODIFIED CLAMP**

**BLEND UP OXYGEN CONSERVER SWITCH AND MONITOR**



Blend Up Kit	**Does Not Include Blender**	BU-017
Flowmeter Bleed On/Off Switch		PRT-2030

**NIF GAUGES / VACUUM PRESSURE GAUGES**

**\*\* -40 to +80 cm H2O GAUGES\*\***

00-274	\$115.00	Manometer with Adjustable Pop Off W/O Max Hand
274-max	\$159.00	Manometer with Pop-Off with MAX Hand
00-310	\$389.00	Manometer with Adjustable Pop-Off 2.5" to 3.5" Clamp
00-300	\$329.00	Manometer with Adjustable Pop-Off 0.75" to 1.75" Clamp
00-199	\$79.00	Replacement Gauge with 1/4" NPT
00-273	\$69.00	Pressure Relief Valve
00-195	\$89.00	Wall Mount Bracket for 00-199



**\*\* ALL 2.5" DIAMETER GAUGES\*\***

**\*\* HEAVY DUTY METAL CASE\*\***



± 30 cm/H2O  
1/8" Tubing Nipple  
30-PBS W/O Recorder  
30-PBR W/Recorder



± 60 cm/H2O  
1/8" Tubing Nipple  
60-PBS W/O Recorder  
60-PBR W/Recorder



± 60 cm/H2O  
1/4" NPT Bottom Port  
60-TBS W/O Recorder  
60-TBR W/Recorder




± 60 cm/H2O  
1/4" NPT Back Port  
60-TRS W/O Recorder  
60-TRR W/Recorder

**GAUGE ADAPTERS**




BE 148-1



BE 148-6



BE 148-4



BE 149-1

**GAUGE COVERS**



\$39.00

#00-150  
Helps Prevent  
Breakage

**MABIS**

**Sphygmomanometer  
Replacement Gauges**



05-230-020 BLACK

05-230-010 BLUE

**GE VERTICAL RAIL  
BRACKET/GAUGE**



**-20 to +100 cm H2O**

\$329.00

GE1072

**± 120 NIF GAUGE**



± 120 cm/H2O  
1/4" NPT Back Port  
120-TRS W/O Recorder  
120-TRR W/Recorder



**Value Med**

**TO PLACE AN ORDER**

**FAX: 888-676-0711**

**EMAIL: valumed.sales@gmail.com**

**SPIROMETERS / CUFFLATORS / PERCUSSORS**



1 1/2"

**Small Wright Mark 8 Respirometer 700-008**



2"

**Large Wright Mark 21 Respirometer 700-021**



**Rubber Crash Guard**

700-401 Small Dial Cover

700-402 Large Dial Face

**CUFF PRESSURE GAUGES**



**STANDARD POSEY CUFFLATOR 0-120 CM H2O 8199**



**POCKET CUFF 0-120 CM H2O 54-04-000**

**Trade In Broken Cufflators & Upgrade for \$179.00**

**SENSITIVE CUFF.....B-54-03-001**

**PULSE OXIMETERS**



**Tuffsat Handheld Pulse Oximeter With 36" Finger Clip Sensor & Carrying Case**

6050-0006-075



**Fingertip Pulse Oximeter Small and Lightweight Bright LED Display**

50-104-004

**36" Tuffsat Finger Clip Replacement**

**Part # 118P46 \$240.000**



**MEDSYSTEMS - Fluid Flow Percussor**

2500 DISS Oxygen w/Hose

2700 DISS Medical Air w/Hose



**G5 NEOCUSSOR Neonatal Percussor Unit**

Neonatal Percussor NC90

Replacement Applicators NAP-7PKG 25/Pack

**TO ORDER: Phone: 888-998-8633 / Fax: 888-676-0711 / Email: Valumed.Sales@Gmail.com**



**REPLACEMENT OXYGEN SENSORS / FUEL CELLS**

**MSA/OHIO**  
406931



12 month warranty  
**Replacement For:**  
MiniOx,  
Esprit Ventilator

**MSA/OHIO**  
472062



12 month warranty  
**Replacement For:**  
Siemens 900C, 300  
Avea Ventilator, PB 840

**MAX-13**



**Replacement For:**  
MiniOX I, II, III, 3000  
Teledyne R-17,  
PB 4-020933-00  
Esprit Ventilator  
14 month warranty

**MAX-13-250**



**Replacement For:**  
MiniOx I, II, III, 3000  
24 month warranty

**MAX-16**



**Replacement For:**  
PB 840 SERIES  
4-072214-00  
16 month warranty

**MAX-12**



**Replacement For:**  
Siemens 900C, 300  
Teledyne R22, R24  
Avea Ventilator  
12 month warranty

**MAX-12A**



**Replacement For:**  
Hillrom Isolette  
Incubator  
Viasys Avea Ventilator  
14 month warranty

**MAX-9**



**Replacement For:**  
Hudson 5500  
12 month warranty

**MAX-1**



**Replacement For:**  
GE/Ohmeda 4760 Oxicap  
5250RGM  
5120 O2 Monitor  
12 month warranty

**MAX-2**



**Replacement For:**  
Draeger anesthesia  
24 month warranty

**MAX-11**



**Replacement For:**  
N.A. Draeger  
Draeger-Narkomed  
Babylog, Evita Ventilator  
12 month warranty

**MAX-12C**



**Replacement For:**  
Viasys, Avea  
Ventilators  
14 month warranty

**MAX-11i**




**Replacement For:**  
Draeger Incubators  
IC 8000  
High Temp Applications  
12 month warranty

**MAX-43**



**Replacement For:**  
GE (Datex Ohmeda)  
Giraffe  
12 month warranty

**MAX-250K**



DC Panduit  
**Replacement For:**  
Viasys, Vela, Ventilators  
DC Adapter.....250K-001  
Panduit.....250K-009

**MAX-17**



**Replacement For:**  
Teledyne  
60T, 191T & 2000T





**Value Med**

**TO PLACE AN ORDER**

**FAX: 888-676-0711**

**EMAIL: valumed.sales@gmail.com**

**OHIO OXYGEN ANALYZERS**

**MAXTEC OXYGEN ANALYZERS**





MiniOX I O2 ANALYZER	
	
473030	\$379.00


MiniOX O2 MONITOR 3000 With HI-LO Alarm	
	
814365	\$603.00


O2 ANALYZER	
	
OM-25AE	\$399.00

O2 MONITOR With HI-LO ALARMS	
	
OM-25ME	\$599.00


**ACCESSORIES**

REPLACEMENT OXYGEN SENSOR		REPLACEMENT SENSOR TEE		REPLACEMENT CABLE			MOUNTING BRACKET			
										
Ohio MiniOx		Maxtec 250E								
OHIO	406931	\$78.00	473021	\$10.50	Ohio	472045	\$59.00	POLE	#474664	\$69.00
MAXTEC	MAX-250E	\$87.00	RP16PO2	\$12.00	Maxtec	R103P16	\$59.00	POLE	R206P75	\$69.00

VENTILATOR SUPPORT ARM WITH FLEXIBLE TIP		
		
116-56	EACH	\$249.00
116-02	2/CASE	\$399.00
116-56K	REPLACEMENT KNOB	\$19.00

SUPPORT ARM BRACKET		
		
Siemens	BE 114-5	\$89.00
Bird/Draeger	RB 002	\$89.00
Vision	RB 002	\$89.00


<b>Replacement For:</b> Vapotherm With O-Ring

SERVO I SUPPORT ARM BRACKET / FISHER & PAYKEL POLE MOUNT BRACKET		
		
9909	EACH	\$69.00
116-99-SE	EACH	\$179.00

**FIBER OPTIC LARYNGOSCOPE BLADES**

MACINTOSH		
PART#	SIZE	PRICE
5-5231-00	Size 0	\$59.00
5-5231-01	Size 1	\$59.00
5-5231-02	Size 2	\$59.00
5-5231-03	Size 3	\$59.00
5-5231-35	Size 3.5	\$59.00
5-5321-04	Size 4	\$59.00

MILLER		
PART#	SIZE	PRICE
5-5233-00	Size 0	\$59.00
5-5233-01	Size 1	\$59.00
5-5233-02	Size 2	\$59.00
5-5233-03	Size 3	\$59.00
5-5233-35	Size 4	\$59.00

WISCONSIN		
PART#	SIZE	PRICE
5-5234-02	Size 2	\$59.00
5-5234-03	Size 3	\$59.00
5-5234-04	Size 4	\$59.00

WIS-HIPPLE		
PART#	SIZE	PRICE
5-5234-20	Size 00	\$59.00
5-5234-00	Size 0	\$59.00
5-5234-01	Size 1	\$59.00
5-5234-15	Size 1.5	\$59.00

FIBER-OPTIC HANDLES		
		<p><b>ALL HANDLES</b> <b>\$68.00</b></p>
PART#	SIZE	
5-0236-09	Medium Handle	"C"
5-0236-10	Penlite Handle	"AA"
5-0236-11	Stubby Handle	"AA"
5-0236-15	Large Handle	"D"
5-0240-52	FO Lamp 5/pack	2.5 V
<p><b>VARIOUS STANDARD AND FIBER OPTIC BLADES TO CHOOSE FROM</b></p>		
*Cranwell	*Phillips	*Robertshaw
*Seward	*Oxford	*Miller/Port
*Heine	*Flangeless Mac	*Siker
*Snow	*Whitehead	*Soper
*Biechman	*Sun Flex	*Prism View

**REPLACEMENT BULBS  
WELCH ALLEN STYLE 5/PACK**



3.5V Laryngoscope  
00300U



3.5V Otoscope  
03100U



2.5V Distal Proctoscope Lamp  
0600U



3.5V Coaxial Ophthalmoscope  
04900U



2.5V Vacuum Otoloscope  
0200U



2.5V Otoscope  
03400U



3.5V Ophthalmoscope Fiber Optic  
03000U



2.5V Ophthalmoscope  
04400U




**Value Med**


**TO PLACE AN ORDER**


**FAX: 888-676-0711**


**EMAIL: [valuemed.sales@gmail.com](mailto:valuemed.sales@gmail.com)**

**STANDARD LARYNGOSCOPE BLADES**

	MACINTOSH		
	PART#	SIZE	PRICE
	5-5052-00	Size 0	\$28.00
	5-5052-01	Size 1	\$28.00
	5-5052-02	Size 2	\$28.00
	5-5052-03	Size 3	\$28.00
	5-5052-04	Size 4	\$28.00

	WISCONSIN		
	PART#	SIZE	PRICE
	5-5057-02	Size 2	\$28.00
	5-5057-03	Size 3	\$28.00
	5-5057-04	Size 4	\$28.00

	MILLER		
	PART#	SIZE	PRICE
	5-5062-20	Size 00	\$28.00
	5-5062-00	Size 0	\$28.00
	5-5062-01	Size 1	\$28.00
	5-5062-15	Size 1.5	\$28.00
	5-5062-02	Size 2	\$28.00
	5-5062-03	Size 3	\$28.00
	5-5062-04	Size 4	\$28.00

	WIS-HIPPLE		
	PART#	SIZE	PRICE
	5-5057-20	Size 00	\$28.00
	5-5057-00	Size 0	\$28.00
	5-5057-01	Size 1	\$28.00
	5-5057-15	Size 1.5	\$28.00

**HANDLES**



PART#	DESCRIPTION		PRICE
5-0237-01	Penlite Handle	“AA”	\$29.50
5-0237-02	Stubby Handle	“AA”	\$29.50
5-0237-03	Medium Handle	“C”	\$29.50
5-0237-05	Large Handle	“D”	\$29.50

**MAGILL FORCEPS**

PART#	DESCRIPTION	PRICE
9-3018-03	Adult	\$16.50
9-3018-02	Child	\$16.50
9-3018-01	Infant	\$16.50



**REPLACEMENT LAMPS - 5/PACK**


*FITS ALL CONVENTIONAL STANDARD BLADES*		
PART#	DESCRIPTION	PRICE
5-0234-01	Standard Frosted, Small	\$17.50
5-0234-04	Standard Frosted, Large	\$17.50
5-0235-01	Standard Clear, Small	\$17.50
5-0235-04	Standard Clear, Large	\$17.50




**SUNBRIGHT REPLACEMENT BULBS**


**\*\*CALL FOR PRICING\*\***


**Welch Allen, Rusch, Propper, Heine Style**

<b>WELCH ALLEN STYLE</b>	
	
00110	Replaced by #04700
00200	2.5V - Vacuum Otoscope
00300	3.5V - Laryngoscope
00400	Replaced by #04800
00600	2.5V - Distal Proctoscope Lamp
00800	6V - Vaginal Speculum
00900	2.5V - Ophthalmoscope - #121
01000	2.5V - Transilluminator
01100	2.5V - Transilluminator
01300	#40000 - Transilluminator
01800	#177 - Reinoscope
02100	2.5V - Ophthalmoscope - #11500
02200	2.5V - Otoscope - #22900
02300	2.5V - Ophthalmoscope - #12900
02500	6V - #46003 - Headlight
02600	6V - #49003 - Headlight
02700	2.5V - Otoscope - #24020
03000	3.5V - Ophthalmoscope - #1160
03000	Type 2 (old style) 3.5V - Ophthalmoscope
03100	3.5V - Otoscope - Fibre Optic
03300	2.5V - Ophthalmoscope - #11511
03400	2.5V - Otoscope - Fibre Optic
03700	3.5V - Retinoscope - #1810
03800	3.5V - Ophthalmoscope - #11810
03900	3.5V - Ophthalmoscope - #11100

<b>WELCH ALLEN STYLE</b>	
	
04000	6V - Sigmioscope
04100	14.25V - Examination Light
04200	12.75V - Examination Light
04300	2.5V - Ophthalmoscope - #11400
04400	2.5V - Ophthalmoscope - Fiber Optic
04500	3.5V - Retinoscope - #18010
04600	2.5V - Retinoscop
04700	2.5V - Laryngoscope Lamp - Small
04800	2.5V - Laryngoscope Lamp - Large
04900	3.5V - Coaxial Ophthalmoscope
06000	2.5V - Fiber Optic Laryngoscope
06000	HPX - Fiber Optic Laryngoscope - 10-32 thread
06100	14.25V - Examination Light
06200	3.5V - Otoscope
06300	20 Watt - Examination Light
06400	12V - Halogen Disply Lamp
06500	3.5V - Otoscope
07200	2.5V - Retinoscope
07800	6V - Vaginal Speculum - #78103
08100	12V - Examination Light
08200	3.5V - Retinoscope - #18200
08300	3.5V - Spot Retinoscope - #18300
08500	3.5V
08800	4.6V

<b>HEINE LAMPS</b>	
034 - 2.5V	051 - 3.5V
035 - 2.5V	058 - 2.5V
037 - 2.5V	059 - 2.5V
038 - 2.5V	060 - 6V
040 - 2.5V	068 - 2.5V
042 - 2.5V	070 - 3.5V
043 - 2.5V	071 - 2.5V
044 - 3.5V	072 - 3.5V
048 - 3.5V	078 - 3.5V
049 - 3.5V	

<b>RUSCH STYLE LAMPS</b>	
	
5630	Lamp for FOCS Laryngoscope
8322	Lamp for Lync fiber optic handle and Optima 2 fibre optic handle
8622 300	Frosted Standard Laryngoscope Lamp - Large
8623 100	Frosted Standard Laryngoscope Lamp - Small

<b>PENLON LAMPS</b>	
	
55091	Laryngoscope - Large
55092	Laryngoscope - Small
55094	Laryngoscope - Fiber Optic



**RESPIRATORY ACCESSORIES**

**TOOTH PROTECTOR**



9-3041-03 Tooth Protector 50/cs

**NOSE CLIPS**



9-0134-44 100/cs

**BACTERIA FILTER/CLEAR**

Bacterial filtration efficiency — 99.8%.  
 Mean particle size 3.1 micron, lab report #18112.  
 22mm Male x 22mm Female, 50/case 64020



**TEST LUNG**

4-0050-50	500-750cc Medium
4-0050-10	1000cc Large



**SILICONE CONNECTOR 22mm FEMALE**

#BE-101-00-3443 SOLD EACH  
 #2200 SOLD IN A CASE OF 50



**MICRO STIM**



8-1053-60 \$239.00

**VENTILATOR ACCESSORIES PB 840**

**PB 840 Ventilator**  
 Disposable Expiratory Bacteria Filter  
 OEM# 4-076887-00 12/Case



**PB 840 Ventilator**  
 Compressor Inlet Filter  
 OEM# 4-07437-00



New PN 10097559  
 Old PN 4-072214-00

**Replacement For:**  
 PB 840  
 Genuine  
 Puritan Bennett  
 Sensor




**7200 Ventilator**

Disposable Expiratory Bacteria Filter  
 OEM# G-060526-00 12/cs



**840 REPLACEMENT BATTERY**




4-070S23-SP

**PB 840 SCREEN LATCH ASSEMBLY**

4-076671-SP DISPLAY		4073020-00 CART	
------------------------	--------------------------------------------------------------------------------------	--------------------	---------------------------------------------------------------------------------------


**BiPAP VISION**



Pollen Filter  
 OEM# 582101 6/Pk.

VISION REPLACEMENT BATTERY 1012819

**BiPAP VISION FILTER COVER**








14-2-23 \$96.00

**V60 VENTILATOR**








Filter 5/pk	V60-6F
Battery	1056921

**ROLLING STANDS**

Standard MRI Compatible	I.V. Stand with 2 or 4 Hook Top	Rolling Stand with Tray	Universal Vent Stand	Notebook Stand
				
WTMRI-100A	1040-2	#1064 12" Tray	#1059	#1268
1059T Std. Pole/Holder	1040-4	#1065 17.5" Tray		

**STAND ACCESSORIES**

Multi-Use Basket	Cylinder Holder	Dual Cylinder Holders		Bag Hook
				
12 x 8 x 6 1034	#1046	#1044	#1044-1	#1047
12 x 12 x 8 1035				

**EMERGENCY PREPAREDNESS**

10 Flowmeter Mounting Bracket With 15' Hose



**HILL ROM BED CYLINDER HOLDER**



1099-BTE

**POLE BRACKET ASSEMBLIES**



Ohio Continuous Suction Regulator With 10' Hose

H-3-PMB

Ohio Intermittent / Continuous Suction Regulator With 10' Hose

H-2-PMB

Oxygen Flowmeter With 10' Hose

15-102-PM



**Value Med**


**TO PLACE AN ORDER**

**FAX: 888-676-0711**

**EMAIL: [valuemed.sales@gmail.com](mailto:valuemed.sales@gmail.com)**


**CYLINDER CARTS & STANDS**

**D & E Size Carts**



PART#	DESCRIPTION	PRICE
B-6105	Standard	\$39.50
B-6105-2U	Standard, Box of 2	\$59.00
B-6C105	Chrome	\$44.00
MRI-1AA	Non-Magnetic	\$179.00
B-6205	Dual Cylinder	\$89.00

**H Cylinder Carts**



SINGLE CART	DESCRIPTION	PRICE
B-6110	Standard	\$289.00
MRI B-6108	Non-Magnetic	\$489.00
DUAL CART	DESCRIPTION	PRICE
B-6210	Standard	\$329.00

**D & E Cylinder Stands**

PART#	DESCRIPTION	PRICE
B-6040	Cylinder Capacity: 4	\$89.00
B-6060	Cylinder Capacity: 6	\$99.00
B-6120	Cylinder Capacity: 12	\$119.00
B-6180	Cylinder Capacity: 18	\$149.00
B-6240	Cylinder Capacity: 24	\$179.00



Multiple D & E Carts & Stands		
2 Wheel		
B-6041	Cylinder Capacity: 4	\$149.00
B-6061	Cylinder Capacity: 6	\$169.00
B-6121	Cylinder Capacity: 12	\$179.50
4 Wheel		
B-6064	Cylinder Capacity: 6	\$196.00
B-6124	Cylinder Capacity: 12	\$219.00
B-6184	Cylinder Capacity: 18	\$229.00
B-6246	Cylinder Capacity: 24	\$269.00





**Disposable Hydrophobic Filter**



**GOMCO**



**AEROS / MOBLVAC / OHIO**

<b>PROTECT SUCTION EQUIPMENT FROM FLUIDS</b>	
01-90-3100	Square 3/Pack
751120	Round W/ Tubing 20/Pack

**OVERFLOW SAFETY TRAP FITTINGS & PARTS**

1. Tubing Nipple	0204-9048-535	
2. Cap	0221-6315-535	
3. O-Ring	0210-0528-300	
4. Splash Deflector	0221-6317-100	
5. Float	0221-6316-100	
6. Jar	0212-0463-100	





**Value Med**

**TO PLACE AN ORDER**

**FAX: 888-676-0711**

**EMAIL: [valuemed.sales@gmail.com](mailto:valuemed.sales@gmail.com)**

**PRECISION MEDICAL SUCTION REGULATORS**

		PRECISION MEDICAL						
		INLET / OUTLET	BARBED HOSE STEM	CHEMETRON	DISS FEMALE	OHMEDA	PURITAN BENNETT	OXEQUIP
2 MODE CONTINUOUS	DISS MALE	PM3021	PM3016	PM3014	PM3015	PM3019	PM3017	PM3018
	HOSE BARB	PM3010	PM3005	PM3003	PM3004	PM3008	PM3006	PM3007
3 MODE CONTINUOUS	DISS MALE	PM3121	PM3116	PM3114	PM3115	PM3119	PM3117	PM3118
	HOSE BARB	PM3110	PM3105	PM3103	PM3104	PM3108	PM3106	PM3107
CONTINUOUS INTERMITTENT	DISS MALE	PM3321	PM3316	PM3314	PM3315	PM3319	PM3317	PM3318
	HOSE BARB	PM3310	PM3305	PM3303	PM3304	PM3308	PM3306	PM3307
PEDIATRIC	DISS MALE	PM3421	PM3416	PM3414	PM3415	PM3419	PM3417	PM3418
	HOSE BARB	PM3410	PM3405	PM3403	PM3404	PM3408	PM3406	PM3407

**\*\*Add PR to Part Number For Preset Precision Medical Version**

**CHEMETRON SUCTION REGULATORS**

		CHEMETRON						
		INLET / OUTLET	DISS FEMALE	BARBED HOSE STEM	OHMEDA	CHEMETRON	OXEQUIP	MEDSTAR
CONTINUOUS/INTERMITTENT	DISS MALE	C-2-130-55	C-2-77-55	C-2-531-55	C-2-431-55	C-2-331-55	C-2-330-55	C-2-231-55
	HOSE BARB	C-2-130-77	C-2-77-77	C-2-531-77	C-2-431-77	C-2-331-77	C-2-330-77	C-2-231-77
CONTINUOUS	DISS MALE	C-3-130-55	C-3-77-55	C-3-531-55	C-3-431-55	C-3-331-55	C-3-330-55	C-3-231-55
	HOSE BARB	C-3-130-77	C-3-77-77	C-3-531-77	C-3-431-77	C-3-331-77	C-3-330-77	C-3-231-77
SURGERY	DISS MALE	C-4-130-55	C-4-77-55	C-4-531-55	C-4-431-55	C-4-331-55	C-4-330-55	C-4-231-55
	HOSE BARB	C-4-130-55	C-4-77-77	C-4-531-77	C-4-431-77	C-4-331-77	C-4-330-77	C-4-231-77
PEDIATRIC CONTINUOUS/INTERMITTENT	DISS MALE	C-5-130-55	C-5-77-55	C-5-531-55	C-5-431-55	C-5-331-55	C-5-330-55	C-5-231-55
	HOSE BARB	C-5-130-55	C-5-77-77	C-5-531-77	C-5-431-77	C-5-331-77	C-5-330-77	C-5-231-77

**\*REMARKETED SUCTION REGULATORS AVAILABLE\***




- \*Chemetron
- \*Ohmeda
- \*Precision Medical
- \*Boehringer
- \*Gentec




**CALL FOR ESTIMATES**

**CANISTERS & ACCESSORIES**

CARDINAL	
	
65651-212	1200cc Stem 40/CS
65651-220	2000cc Stem 40/CS
65651-230	3000cc Stem 40/CS

BEMIS	
	
424410	800cc Stem 100/CS
434410	1200cc Stem 48/CS
435410	2000cc Stem 48/CS

ALLIED	
	
08-003	1500ML Stem 48/Case
08-004	1500ML DISS 48/Case
08-002	2400ML Stem 36/Case
08-001	2400ML DISS 36/Case

**REUSABLE OUTER CANISTER FOR DISPOSABLE LINERS**

*BUILT IN* ON/OFF VALVE	
	
65652-511	1000cc
65652-516	1500cc
65652-531	3000cc

WALL MOUNT	
	
65652-513	1000cc
65652-518	1500cc
65652-533	3000cc

EXTENDED BRACKET FOR REGULATORS	
	
65652-512	1000cc
65652-517	1500cc
65652-532	3000cc

DIRECT TO REGULATOR	
	
65652-555	1000cc
65652-562	1500cc

**CARDINAL**

Canister Ring Featuring Pre-Attached Tubing & Built-In Internal On/Off Valve  
 Accommodates 1200cc Canister with On/Off Valve **65652-159**  
 Accommodates 2000cc & 3000cc Canister with On/Off Valve **65652-139**  
 Accommodates 1200cc Canister without On/Off Valve **65652-155**  
 Accommodates 2000cc & 3000cc Canister without On/Off Valve **65652-135**



**DISPOSABLE SUCTION TUBING**


Suction Tubing, 6 Ft 536010 25/Cs  
 Suction Tubing, 1.5 Ft 536510 50/Cs



**SUCTION CANISTER RINGS**

**BEMIS WALL CANISTER RING - METAL**

Accommodates Bemis 800cc Canister	4 3/4"ID	534510
Accommodates Bemis 1200, 2000 & 3000cc Canisters	5 5/8"ID	532810
Accommodates Medi-Vac 1200cc Guardian Canister	4 3/4"ID	534510



**BEMIS WALL CANISTER RING - PLASTIC**

Accommodates 1200cc, 2000cc & 3000cc Bemis Canister	5 5/8"ID	532710
-----------------------------------------------------	----------	--------



**ALLIED HEALTHCARE CANISTER RING**

Accommodates 1500 & 2400ml Allied Canisters	5 1/4"ID	02-0150
---------------------------------------------	----------	---------






**CARDINAL STANDARD CANISTER RING**

-Accommodates Cardinal 1000cc Canister	4 3/8"ID	65652-584
-Accommodates Cardinal 1200cc Canister	4 3/8"ID	65652-584
-Accommodates Cardinal 1500cc Canister	5 1/8"ID	65652-563
-Accommodates Cardinal 3000cc Canister	6 3/4"ID	65652-561
-Accommodates Medi-Vac 1200cc	4 3/8"ID	65652-584
& Guardian / Bemis 800cc Canister	4 3/8"ID	65652-584





**FOUR CANISTER CARRIERS**


**\*\*STANDS COME WITHOUT RINGS\*\***

ALLIED / CHEMETRON			BEMIS			CARDINAL		
								
20-02-0023	36" Tall	\$313.00	533910	ADJUSTABLE	\$313.00	652-586	30" Tall	\$313.00

**OHIO / AMVEX / GENTEC SUCTION REGULATORS**

		OHIO / OHMEDA PUSH TO SET								
		INLET OUTLET	DISS FEMALE	BARBED HOSE STEM	OHMEDA	CHEMETRON	OXEQUIP	MEDSTAR	PURITAN BENNETT	
	CONTINUOUS/ INTERMIT- TENT	DISS MALE	HP-2-130-55	HP-2-77-55	HP-2-531-55	HP-2-431-55	HP-2-331-55	HP-2-330-55	HP-2-231-55	
		HOSE BARB	HP-2-130-77	HP-2-77-77	HP-2-531-77	HP-2-431-77	HP-2-331-77	HP-2-330-77	HP-2-231-77	
	CONTINUOUS	DISS MALE	HP-3-130-55	HP-3-77-55	HP-3-531-55	HP-3-431-55	HP-3-331-55	HP-3-330-55	HP-3-231-55	
		HOSE BARB	HP-3-130-77	HP-3-77-77	HP-3-531-77	HP-3-431-77	HP-3-331-77	HP-3-330-77	HP-3-231-77	
	SURGERY	DISS MALE	HP-4-130-55	HP-4-77-55	HP-4-531-55	HP-4-431-55	HP-4-331-55	HP-4-330-55	HP-4-231-55	
		HOSE BARB	HP-4-130-55	HP-4-77-77	HP-4-531-77	HP-4-431-77	HP-4-331-77	HP-4-330-77	HP-4-231-77	
<b>PUSH TO SET CONT / INT</b>	PEDIATRIC CONTINUOUS/ INTERMIT- TENT	DISS MALE	HP-5-130-55	HP-5-77-55	HP-5-531-55	HP-5-431-55	HP-5-331-55	HP-5-330-55	HP-5-231-55	
		HOSE BARB	HP-5-130-55	HP-5-77-77	HP-5-531-77	HP-5-431-77	HP-5-331-77	HP-5-330-77	HP-5-231-77	

		AMVEX								
		INLET OUTLET	DISS FEMALE	BARBED HOSE STEM	OHMEDA	CHEMETRON	OXEQUIP	MEDSTAR	PURITAN BENNETT	
	CONTINUOUS/ INTERMITTENT	DISS MALE	A-2-130-55	A-2-77-55	A-2-531-55	A-2-431-55	A-2-331-55	A-2-330-55	A-2-231-55	
		HOSE BARB	A-2-130-77	A-2-77-77	A-2-531-77	A-2-431-77	A-2-331-77	A-2-330-77	A-2-231-77	
	CONTINUOUS	DISS MALE	A-3-130-55	A-3-77-55	A-3-531-55	A-3-431-55	A-3-331-55	A-3-330-55	A-3-231-55	
		HOSE BARB	A-3-130-77	A-3-77-77	A-3-531-77	A-3-431-77	A-3-331-77	A-3-330-77	A-3-231-77	
	SURGERY	DISS MALE	A-4-130-55	A-4-77-55	A-4-531-55	A-4-431-55	A-4-331-55	A-4-330-55	A-4-231-55	
		HOSE BARB	A-4-130-55	A-4-77-77	A-4-531-77	A-4-431-77	A-4-331-77	A-4-330-77	A-4-231-77	
<b>AMVEX CONT / INT</b>	PEDIATRIC CONTINUOUS/ INTERMITTENT	DISS MALE	A-5-130-55	A-5-77-55	A-5-531-55	A-5-431-55	A-5-331-55	A-5-330-55	A-5-231-55	
		HOSE BARB	A-5-130-55	A-5-77-77	A-5-531-77	A-5-431-77	A-5-331-77	A-5-330-77	A-5-231-77	

		GENTEC								
		INLET OUTLET	DISS FEMALE	BARBED HOSE STEM	OHMEDA	CHEMETRON	OXEQUIP	MEDSTAR	PURITAN BENNETT	
	CONTINUOUS/ INTERMITTENT	DISS MALE	G-2-130-55	G-2-77-55	G-2-531-55	G-2-431-55	G-2-331-55	G-2-330-55	G-2-231-55	
		HOSE BARB	G-2-130-77	G-2-77-77	G-2-531-77	G-2-431-77	G-2-331-77	G-2-330-77	G-2-231-77	
	CONTINUOUS	DISS MALE	G-3-130-55	G-3-77-55	G-3-531-55	G-3-431-55	G-3-331-55	G-3-330-55	G-3-231-55	
		HOSE BARB	G-3-130-77	G-3-77-77	G-3-531-77	G-3-431-77	G-3-331-77	G-3-330-77	G-3-231-77	
	SURGERY	DISS MALE	G-4-130-55	G-4-77-55	G-4-531-55	G-4-431-55	G-4-331-55	G-4-330-55	G-4-231-55	
		HOSE BARB	G-4-130-55	G-4-77-77	G-4-531-77	G-4-431-77	G-4-331-77	G-4-330-77	G-4-231-77	
<b>GENTEC CONT / INT</b>	PEDIATRIC CON- TINUOUS/ INTERMITTENT	DISS MALE	G-5-130-55	G-5-77-55	G-5-531-55	G-5-431-55	G-5-331-55	G-5-330-55	G-5-231-55	
		HOSE BARB	G-5-130-55	G-5-77-77	G-5-531-77	G-5-431-77	G-5-331-77	G-5-330-77	G-5-231-77	






**Value Med**

**TO PLACE AN ORDER**

**FAX: 888-676-0711**




**EMAIL: [valuemed.sales@gmail.com](mailto:valuemed.sales@gmail.com)**

**OHIO / OHMEDA / SUCTION REGULATORS**

	OHIO / OHMEDA								
	MODEL	INLET / OUTLET	DISS FEMALE	BARBED HOSE STEM	OHMEDA	CHEMETRON	OXEQUIP	MEDSTAR	PURITAN BENNETT
CONTINUOUS/ INTERMITTENT	DISS MALE	H-2-130-55	H-2-77-55	H-2-531-55	H-2-431-55	H-2-331-55	H-2-330-55	H-2-231-55	
	HOSE BARB	H-2-130-77	H-2-77-77	H-2-531-77	H-2-431-77	H-2-331-77	H-2-330-77	H-2-231-77	
CONTINUOUS	DISS MALE	H-3-130-55	H-3-77-55	H-3-531-55	H-3-431-55	H-3-331-55	H-3-330-55	H-3-231-55	
	HOSE BARB	H-3-130-77	H-3-77-77	H-3-531-77	H-3-431-77	H-3-331-77	H-3-330-77	H-3-231-77	
SURGERY	DISS MALE	H-4-130-55	H-4-77-55	H-4-531-55	H-4-431-55	H-4-331-55	H-4-330-55	H-4-231-55	
	HOSE BARB	H-4-130-55	H-4-77-77	H-4-531-77	H-4-431-77	H-4-331-77	H-4-330-77	H-4-231-77	
OHIO / OHMEDA CONT / INT	PEDIATRIC CONTINUOUS/ INTERMITTENT	DISS MALE	H-5-130-55	H-5-77-55	H-5-531-55	H-5-431-55	H-5-331-55	H-5-330-55	H-5-231-55
	HOSE BARB	H-5-130-55	H-5-77-77	H-5-531-77	H-5-431-77	H-5-331-77	H-5-330-77	H-5-231-77	

**\*\*\*\*\*Protect Your Suction Equipment\*\*\*\*\***

**SUCTION REGULATOR OVERFLOW SAFETY TRAP**

<p>DISS Female</p> 	<p>ADD "T" TO SUCTION REGULATOR PART NUMBER FOR TRAP</p> 	<p>Locking Gland</p> 
<p>130-1000      \$63.50</p>		<p>295-1000      \$69.50</p>

**Overflow Safety Trap Features And Benefits**

- ✓ Limit the contamination risk of the regulator and the vacuum pipeline network.
- ✓ Unique float design ensures maximum vacuum flow and positive shutoff.
- ✓ Clear, shatter resistant polycarbonate bottle can be steam or sterilized after using.
- ✓ It helps prevent aspirate from entering the regulator as a result of misuse and extends the life of suction equipment and pipeline systems.
- ✓ Meets NFPA and ASTM standards.

**TO ORDER: Phone: 888-998-8633 / Fax: 888-676-0711 / Email: [Valuemed.Sales@Gmail.com](mailto:Valuemed.Sales@Gmail.com)**

**SUCTION PUMPS**

**BATTERY OPERATED CRASH CART PUMPS**



MANUFACTURER	PART#	RUN TIME
GOMCO	G-180	1 Hour
PRECISION MEDICAL	PM 65	75 Minutes
OHIO	758000	1 Hour

**TABLE TOP SUCTION PUMPS**



MANUFACTURER	MODEL	PART#
GOMCO	Model 300	01-22-3001
GOMCO	Model 405	01-12-4005

**STAND MODEL SUCTION PUMPS**



MANUFACTURER	MODEL	PART#
GOMCO	3040 Aspirator	01-32-3041
GOMCO	6000 Thermotic	01-32-6002
GOMCO	6037 Combination Unit	01-32-6037
AEROS	Mobilvac III Combination Unit	70-850-01



**\*SUCTION REGULATOR REPAIR\***

\* Chemetron \* Ohmeda \* Precision Medical \* Boehringer \* Gentec\*

**ALL MODELS AVAILABLE: (CONTINUOUS-INTERMITTENT-SURGERY)**






**With Suction Regulator  
Repair or Upgrade**






**Your Overflow Safety Trap is  
\$39.50**

**SUCTION ACCESSORIES**

<b>CEILING MOUNT 7" LONG ADAPTER</b>



<b>**Adapter Only** 852</b>

	<b>1/8" NPT</b>
	DISS Male, Vacuum
25-132-02	\$17.00
	<b>Tubing Nipple</b>
	Vacuum Regulator
25-999	\$9.50
	<b>Disposable</b>
	Hydrophobic Filter 3/Pack
01-90-3100	\$19.00


	<b>Metal Female DISS x Universal Tubing Barb</b>
531410	\$29.00
	<b>Plastic Female DISS with 1/8" NPT Adapter</b>
65652-131	\$18.50
	<b>Female DISS x Tubing Barb</b>
531110	\$18.75
	<b>1/8" Male NPT x DISS x Universal Tubing Barb</b>
531110	\$27.50
	<b>Main Line Suction Full Demand</b>
	<b>INLET      PART#</b>
	DISS Female Inlet      40-131
	Ohmeda Inlet      40-531
	Chemetron / NCG Inlet      40-431
	Oxequip Inlet      40-331
	Medstar Inlet      40-332
	Puritan Bennett Inlet      40-231

<b>OVERFLOW FILTER SET-UP</b>



3100-01 Flexi-Filter
3100-44 Flexi-Filter w/Elbow
\$37.50 Each


<b>SUCTION REGULATOR SECURITY CABLE</b>	
	
#OH-167	\$5.35 Each

**BRACKETS**

<b>VERTICAL RAIL BRACKET</b>

\$259.00 EACH
Bed Warmers / OR Stands
OXYGEN      #6700-0411-800
MED AIR      #6700-0412-800
VACUUM      #6700-0410-800

<b>INCLINE VACUUM ON/OFF VALVE</b>	
	Plastic
#45-65652-580	\$18.50
	Metal
1/8" NPT	
#45-7064	\$49.00


<b>Fairfield / HilRom Rail Mount</b>	<b>Bottle Holder</b>
	
#40-9788      \$129.00	#22-11-0004      \$39.00

<b>RAIL BLOCK ASSEMBLY</b>

OXYGEN      7223
MED AIR      7224
VACUUM      7225













**WALL MOUNT SUCTION CANISTER HOLDERS**

	Plastic		Aluminum
530510	\$19.50	#22-11-0006	\$49.50

**VACUUM QUICK CONNECT WITH BOTTLE HOLDER**

		
<b>INLET</b>	<b>PART#</b>	<b>PRICE</b>
DISS Female Inlet	40-30-131	\$139.00
Ohmeda Inlet	40-30-531	\$139.00
Chemetron / NCG Inlet	40-30-431	\$139.00
Oxequip Inlet	40-30-331	\$139.00
Medstar Inlet	40-30-332	\$159.00
Puritan Bennett Inlet	40-30-231	\$139.00

**SPECIALTY FLOWMETERS**

MODEL						
	O2 NEB BLOCK 0-15 FLOWMETER 7-8,14-16,24-26	O2 ELIMINATOR 0-15LPM FLOWMETER	O2 SELECT FLOWMETER 0-15 LPM	O2 DIAL FLOWMETER 25,15,8 & 4 LPM	AIR DIAL FLOWMETER 25 & 15 LPM	CO2 0-12LPM
DISS FEMALE 	1MFA6073 7-8LPM 1MFA6153 14-16LPM 1MFA6243 24-26LPM	1MFA7153 15LPM 1MFA7243 24LPM 1MFA7002 FLUSH	3MFA-1003	25-100-D 25 LPM 15-100-D 15LPM 8-100-D 8LPM 4-100-D 4LPM	25-110-D 25 LPM 15-110-D 15LPM 8-110-D 8LPM 4-110-D 4LPM	CO2-180
PURITAN BENNETT 	1MFA6078 7-8LPM 1MFA6158 14-16LPM 1MFA6248 24-26LPM	1MFA7158 15LPM 1MFA7248 24LPM 1MFA7008 FLUSH	3MFA-1008	25-201-D 25 LPM 15-201-D 15LPM 8-201-D 8LPM 4-201-D 4LPM	25-210-D 25 LPM 15-210-D 15LPM 8-210-D 8LPM 4-210-D 4LPM	CO2-280
CHEMETRON 	1MFA6076 7-8LPM 1MFA6156 14-16LPM 1MFA6246 24-26LPM	1MFA7156 15LPM 1MFA7246 24LPM 1MFA7006 FLUSH	3MFA-1006	25-401-D 25 LPM 15-401-D 15LPM 8-401-D 8LPM 4-401-D 4LPM	25-410-D 25 LPM 15-410-D 15LPM 8-410-D 8LPM 4-410-D 4LPM	CO2-480
OXEQUIP 	1MFA6077 7-8LPM 1MFA6157 14-16LPM 1MFA6247 24-26LPM	1MFA7157 15LPM 1MFA7247 24LPM 1MFA7007 FLUSH	3MFA-1007	25-301-D 25 LPM 15-301-D 15LPM 8-301-D 8LPM 4-301-D 4LPM	25-311-D 25 LPM 15-311-D 15LPM 8-311-D 8LPM 4-311-D 4LPM	N/A
MEDSTAR 	1MFA607M 7-8LPM 1MFA615M 14-16LPM 1MFA624M 24-26LPM	1MFA715M 15 LPM 1MFA724M 24LPM 1MFA700M FLUSH	3MFA-100M	25-302-D 25 LPM 15-302-D 15LPM 8-302-D 8LPM 4-302-D 4LPM	25-312-D 25 LPM 15-312-D 15LPM 8-312-D 8LPM 4-312-D 4LPM	CO2-355
OHMEDA 	1MFA6075 7-8LPM 1MFA6155 14-16 LPM 1MFA6245 24-26LPM	1MFA7155 15LPM 1MFA7245 24LPM 1MFA7005 FLUSH	3MFA-1005	25-501-D 25 LPM 15-501-D 15LPM 8-501-D 8LPM 4-501-D 4LPM	25-510-D 25 LPM 15-510-D 15LPM 8-510-D 8LPM 4-510-D 4LPM	CO2-580
1/8" NPT	1MFA6070 7-8LPM 1MFA6150 14-16LPM 1MFA6240 24-26LPM	1MFA7150 15LPM 1MFA7240 24LPM 1MFA7000 FLUSH	3MFA-1001	18-25-22D 25LPM 18-15-22D 15LPM 18-8-22D 8LPM 18-4-22D 4LPM	18-25-23D 25LPM 18-15-23D 15LPM 18-8-23D 8LPM 18-4-23D 4LPM	18-CO2

















**Value Med**

**TO PLACE AN ORDER**


**FAX: 888-676-0711**


**EMAIL: [valuemed.sales@gmail.com](mailto:valuemed.sales@gmail.com)**

## SPECIALTY FLOWMETERS


MODEL						
DISS FEMALE 	TM-100	NF-100	CM-100	N20-120	75-100	75-110
PURITAN BENNETT 	TM-201	NF-201	CM-201	N20-220	75-201	75-210
CHEMETRON 	TM-401	NF-401	CM-401	N20-420	75-401	75-410
OXEQUIP 	TM-301	NF-301	CM-301	N20-321	75-301	75-311
MEDSTAR 	TM-302	NF-302	CM-302	N20-322	75-302	75-312
OHMEDA 	TM-501	NF-501	CM-501	N20-520	75-501	75-510
1/8" NPT	18-TM	18-NF	18-CM	18-N20	18-75-22	18-75-23

**OXYGEN \* MEDICAL AIR \* NITROUS OXIDE  
HIGH PRESSURE REGULATORS**

50 PSI PRESET REGULATORS		
		
<b>OXYGEN</b>		<b>\$159.00</b>
<b>50-870</b>	Small Cylinder E Yoke	870 CGA
<b>50-540</b>	Large Cylinder Nut	540 CGA
<b>AIR</b>		<b>\$179.00</b>
<b>50-950-08</b>	Small Cylinder E Yoke	950 CGA
<b>50-346-08</b>	Large Cylinder Nut	346 CGA
<b>NITROUS OXIDE</b>		<b>\$289.00</b>
<b>50-910</b>	Small Cylinder E Yoke	910 CGA
<b>50-326</b>	Large Cylinder Nut	326 CGA


0-100 PSI ADJUSTABLE REGULATORS		
		
<b>OXYGEN</b>		<b>\$199.00</b>
<b>100-870</b>	Small Cylinder E Yoke	870 CGA
<b>100-540</b>	Large Cylinder Nut	540 CGA
<b>AIR</b>		<b>\$249.00</b>
<b>100-950</b>	Small Cylinder E Yoke	950 CGA
<b>100-346</b>	Large Cylinder Nut	346 CGA
<b>NITROUS OXIDE</b>		<b>\$289.00</b>
<b>100-910</b>	Small Cylinder E Yoke	910 CGA
<b>100-326</b>	Large Cylinder Nut	326 CGA

**HELI-OX HIGH PRESSURE REGULATORS \$329.00**


0-100 PSI ADJUSTABLE		
		
<b>100% HELIUM</b>		
<b>100-930</b>	Small Cylinder E Yoke	930 CGA \$329.00
<b>100-580</b>	Large Cylinder Nut	580 CGA \$329.00
<b>HELIUM &lt;80%</b>		
<b>100-890</b>	Small Cylinder E Yoke	880 CGA \$329.00
<b>100-280</b>	Large Cylinder Nut	280 CGA \$329.00
<b>HELIUM &gt; 80%</b>		
<b>100-930</b>	Small Cylinder E Yoke	930 CGA \$329.00
<b>100-500-1180</b>	Lrge Cylinder Nut	500 CGA \$329.00




**REPLACEMENT**  
0-16 LPM  
HELI-OX  
FLOWMETER  
**Part # 65-HE**  
**\$139.00**

50 PSI PRESET		
		
<b>100% HELIUM</b>		
<b>50-930-65</b>	Small Cylinder E Yoke	930 CGA \$329.00
<b>50-580-65</b>	Large Cylinder Nut	580 CGA \$329.00
<b>HELIUM &lt;80%</b>		
<b>50-890-65</b>	Small Cylinder E Yoke	880 CGA \$329.00
<b>50-280-65</b>	Large Cylinder Nut	280 CGA \$329.00
<b>HELIUM &gt; 80%</b>		
<b>50-930-65</b>	Small Cylinder E Yoke	930 CGA \$329.00
<b>50-500-65</b>	Large Cylinder Nut	500 CGA \$329.00




**CARBON DIOXIDE HIGH PRESSURE REGULATORS**

CARBON DIOXIDE REGULATORS WITH FLOWMETERS			
<b>100% CARBON DIOXIDE</b>			Replacement CO2 Flowmeter 0-12LPM  <b>\$139.00</b> Starts at 1/2 Liter Part # 18-CO2
0-12 LPM Flowmeter			
50-940-18 Small Cylinder E Yoke	940 CGA	\$329.00	
50-320-18 Large Cylinder Nut	320 CGA	\$329.00	
<b>CARBON DIOXIDE &gt; 7% &amp; OXYGEN</b>			
0-12 LPM Flowmeter			
50-940-18 Small Cylinder E Yoke	940 CGA	\$329.00	
50-500-18 Large Cylinder Nut	500 CGA	\$329.00	
<b>CARBON DIOXIDE &lt; 7% &amp; OXYGEN</b>			
0-12 LPM Flowmeter			
50-880-18 Small Cylinder E Yoke	880 CGA	\$329.00	
50-280-18 Large Cylinder Nut	380 CGA	\$329.00	

50 PSI PRESET CARBON DIOXIDE REGULATORS		
		
<b>100% CARBON DIOXIDE</b>		
50-940 Small Cylinder E Yoke	940 CGA	\$329.00
50-320 Large Cylinder Nut	320 CGA	\$329.00
<b>CARBON DIOXIDE &lt;7% &amp; OXYGEN</b>		
50-880 Small Cylinder E Yoke	880 CGA	\$329.00
50-280 Large Cylinder Nut	380 CGA	\$329.00
<b>CARBON DIOXIDE &gt; 7% &amp; OXYGEN</b>		
50-940 Small Cylinder E Yoke	940 CGA	\$329.00
50-500-1200 Large Cylinder Nut	500 CGA	\$329.00

ADJUSTABLE CARBON DIOXIDE REGULATORS					
ADJUSTABLE 0-100 PSI		ADJUSTABLE 0-30 PSI With Tubing Nipple Clinical Instrumentation			
					
<b>100% CARBON DIOXIDE</b>		<b>100% CARBON DIOXIDE</b>			
100-940 Small Cylinder E Yoke	940 CGA	\$329.00	30-940 Small Cylinder E Yoke	940 CGA	\$329.00
100-320 Large Cylinder Nut	320 CGA	\$329.00	30-320 Large Cylinder Nut	320 CGA	\$329.00
<b>CARBON DIOXIDE &lt; 7% &amp; OXYGEN</b>		<b>CARBON DIOXIDE &lt;7% &amp; OXYGEN</b>			
100-880 Small Cylinder E Yoke	880 CGA	\$329.00	30-880 Small Cylinder E Yoke	880 CGA	\$329.00
100-280 Large Cylinder Nut	380 CGA	\$329.00	30-280 Large Cylinder Nut	380 CGA	\$329.00
<b>CARBON DIOXIDE &gt; 7% &amp; OXYGEN</b>		<b>CARBON DIOXIDE &gt; 7% &amp; OXYGEN</b>			
100-940 Small Cylinder E Yoke	940 CGA	\$329.00	30-940 Small Cylinder E Yoke	940 CGA	\$329.00
100-500-1200 Large Cylinder Nut	500 CGA	\$329.00	30-500-1200 Large Cylinder Nut	500 CGA	\$329.00

**HIGH PRESSURE NITROGEN REGULATORS**

								
Adjustable, 0-200 PSI, Diaphragm			Adjustable, 0-200 PSI, Diaphragm			Adjustable, 0-200 PSI, Diaphragm		
2" Gauges with DISS Outlet			2" Gauges & Schrader Female			2" Gauge, Schrader Coupler & DISS		
\$389.00			\$389.00			\$489.00		
200-960-142	Small Cylinder E Yoke	960 CGA	200-960-610	Small Cylinder E Yoke	960 CGA	200-960-610-142	Small E Cylinder Yoke	960 CGA
200-580-142	Large Cylinder Nut	580 CGA	200-580-610	Large Cylinder Nut	580 CGA	200-580-610-142	Large Cylinder Nut	580 CGA

**LOW PRESSURE NITROGEN REGULATORS \$469.00**

WALL MOUNTED	
	
Elbowed Schrader Male to Schrader Female Coupler	615-610
Elbowed Schrader Male to DISS Male Coupler	615-142
Elbowed DISS Female to DISS Male	144-142
Elbowed DISS Female to Schrader Female	144-610

COLUMN MOUNTED	
	
Schrader Male to Schrader Female Coupler	612-610
Schrader Male to DISS Male Coupler	612-142
DISS Female to DISS Male	140-142
DISS Female to Schrader Female	140-610


N2 Regulator With Flowmeter High Pressure	
Large Cylinder	50-580-112
Small Cylinder	50-960-112

REPLACEMENT REGULATOR KNOB ASSEMBLY	
	
Assembly	\$239.00 11-91-0170
Knob Only	\$69.00 11-91-0172

Davol Arthro-Flo Regulator	
	
0037760	\$925.00







**NITROGEN ACCESSORIES**


	NITROGEN HOSE ASSEMBLIES \$239.00		
	611-20-26-610	20 FT Nitrogen Hose	Schrader Male x Schrader Female
	140-20-26-610	20 FT Nitrogen Hose	DISS Nitrogen x Schrader Female
	140-20-26-140	20 FT Nitrogen Hose	DISS Nitrogen x DISS Nitrogen



	NITROGEN FLEXIBLE CONVERTER HOSE	
	Air for Nitrogen to DISS Quick Connect	
	55-1-610	

<b>Schrader Elbow Quick Connect</b>	
55-E-610	


	NITROGEN WYE ASSEMBLIES \$289.00		
	DISS Female x DISS Male	140-142-142	
	DISS Female x Schrader Female	140-610-610	
	Schrader Male x Schrader Female	612-610-610	
	Schrader Male x DISS Male	612-142-142	

SCHRADER MALE QUICK CONNECT	
Air / Nitrogen 1/8" NPT	
	
Swivel Style Quick Connect N2	611-18
Non Swivel Quick Connect N2	612-18

DISS Female Nitrogen 1/4" NPT		
Hand Tight	140-14	
Hex Nut	141-14	

	SCHRADER FEMALE COUPLER AIR FOR NITROGEN 1/4" NPT	
	610-14	
	DISS MALE NITROGEN 1/4" NPT	
	DISS Male with Check	142-14-CV
	DISS Male without Check	142-14

BLACK NITROGEN 2 PIECE KIT	
For Nitrogen Only	
	
640-641-2	

	BLACK SCHRADER NITROGEN WYE	
	DISS Inlet	22-640-140
	Schrader Inlet	22-641-641

**INSUFFLATOR / LUNG DIFFUSION REGULATOR**

**DISS x DISS Female**  
5Ft CO2 200 PSI



180-5-28-180

**Insufflator Yoke**




94-06

**Male Connector**




19-276

**Insufflator High Pressure Hose**  
3000 PSI




716-3-716

**High Pressure Hose 2 Tank Hook-Up 3000 psi Max With 2 Contents Gauges**



G2-716-31-03

**E TANK YOKE**




PART#	GAS	CGA
870-14	Oxygen	870
880-14	CO2-O2	880
910-14	N2O	910
940-14	CO2	940
950-14	AIR	950
960-14	N2	960
973-14	GAS MIX	973
930-14	HELI-OX	930
890-14	GAS MIX	890

**H-TANK HEX NUT AND NIPPLE**



PART#	GAS	CGA
280-14	HELI-OX	280
236-14	N2O	236
346-14	AIR	346
500-14	GAS MIX	500
540-14	OXYGEN	540
580-14	N2 & HELI	580
320-14	CO2	320


**LUNG DIFFUSION REGULATOR**




<b>30-500-LD</b>	30 PSI W/ NIPPLE	500 CGA	\$369.00
<b>60-500-LD</b>	60 PSI W/ NIPPLE	500 CGA	\$369.00


**FLEXIBLE STAINLESS STEEL PIGTAILS**

\*OXYGEN \* MEDICAL AIR \* NITROUS OXIDE \* NITROGEN \* CARBON DIOXIDE \*


 \$239.00 3000 PSIG	FOR MANIFOLD REPLACEMENT		Flexible Stainless Steel Pigtail
	326-32-326	CGA 326	Nitrous Oxide
	346-32-346	CGA 346	Medical Air
	580-32-580	CGA 580	Nitrogen
	540-31-540	CGA 540	Oxygen (Copper Only)
	320-32-320	CGA 320	Carbon Dioxide

 Custom Configurations Available Upon Request 3000 PSIG \$349.00	MULTIPLE CYLINDER TANK HOOKUP		
			30" Long with Coupler Tee
	T326-30-326	CGA 326	Nitrous Oxide
	T346-30-346	CGA 346	Medical Air
	T580-30-580	CGA 580	Nitrogen
	T540-31-540	CGA 540	Oxygen (Copper Only)
	T320-30-320	CGA 320	Carbon Dioxide

**WALL MOUNT RAIL SYSTEM**

* FAIRFIELD * HOSPITAL SYSTEMS * EASTERN * HILL-ROM *	
	
*CHEMETRON	*OXYGEN
*OHMEDA	*MEDICAL AIR
*PURITAN BENNETT	*VACUUM
*MEDSTAR	*NITROUS OXIDE
*DISS	
DISS Standard 2 Foot Hose - Custom Lengths Available	
R1 for 1 Outlet - R2 for 2 Outlets / Flowmeters	

**WALL MOUNT EXTENSION MANIFOLD**

	
*CHEMETRON	*OXYGEN
*OHMEDA	*MEDICAL AIR
*PURITAN BENNETT	*VACUUM
*MEDSTAR	*NITROUS OXIDE
*DISS	
Standard 2 Foot Hose - Custom Lengths Available	

**HORIZON HEADWALL SYSTEMS MODULE**

HOSE FROM THE BOTTOM	
CHEMETRON	
Oxygen	Z-42-400
Medical Air	Z-42-410
Vacuum	Z-42-430
OHIO	
Oxygen	Z-42-500
Medical Air	Z-42-510
Vacuum	Z-42-530
PURITAN BENNETT	
Oxygen	Z-42-200
Medical Air	Z-42-210
Vacuum	Z-42-230
DISS	
Oxygen	Z-42-102
Medical Air	Z-42-112
Vacuum	Z-42-132

HOSE FROM THE BACK	
CHEMETRON	
Oxygen	ZB-42-400
Medical Air	ZB-42-410
Vacuum	ZB-42-430
OHIO	
Oxygen	ZB-42-500
Medical Air	ZB-42-510
Vacuum	ZB-42-530
PURITAN BENNETT	
Oxygen	ZB-42-200
Medical Air	ZB-42-210
Vacuum	ZB-42-230
DISS	
Oxygen	ZB-42-102
Medical Air	ZB-42-112
Vacuum	ZB-42-132

**HORIZON REPLACEMENT CLIPS**

**\*\*\*INCLUDES 3 FOOT DISS CONNECTION HOSE\*\*\***



Oxygen #877 CLIP ONLY	Medical Air #878 CLIP ONLY	Vacuum #876 CLIP ONLY
\$44.50 Each	\$44.50 Each	\$44.50 Each

**ENVIROVAC**



64-90-VAC

Outlet Cleaner


**ENVIROVAC**

- ✓ ENVIROVAC is an enzyme based product designed to remove debris from vacuum system outlets and all subsidiary vacuum components.
- ✓ This process can eliminate the need to isolate, disassemble and rebuild outlets.
- ✓ Increases restricted vacuum flow.
- ✓ Easy and safe to use.
- ✓ Safe for all materials: hoses, piping, valves and O-rings.

**HIGH PRESSURE HOSES AND HOSE ACCESSORIES**






HALO TUBING SUPPORT	
PART#	SIZE
4310-02	Size 2 W/ Std. Base
4310-03	Size 3 W/ Std. Base
4310-04	Size 4 W/ Std. Base

AIR HOSE CRIMPER	CR-100
	

AIR GUN WITH COILED HOSE - 20 Ft


OUTLET CONNECTION		
Ohio Medical Air	511-40	\$139.00
Chemetron Medical Air	411-40	\$139.00
Puritan Bennett Medical Air	211-40	\$139.00
DISS Medical Air	112-40	\$139.00
DISS Instrument Air	277-40	\$279.00
DISS Nitrogen	140-40	\$279.00

COILED HIGH PRESSURE HOSES		
		
DUAL COILED HOSE OXYGEN/AIR		
Dual 10' Blender	101-D10-111	\$195.00
Dual 15' Ventilator	100-D15-110	\$229.00
		
SINGLE HOSE 15 FOOT LENGTH		
Oxygen	100-015-100	\$119.00
Air	110-A15-110	\$119.00

HIGH PRESSURE AIR GUN WITH REGULATOR - 20 Ft		
		
Air Gun with Coiled Hose and Adjustable Regulator for Wall Outlets		
OUTLET CONNECTION		
Ohio Medical Air	511-42	\$249.00
Chemetron Medical Air	411-42	\$249.00
Puritan Bennett Medical Air	211-42	\$249.00
DISS Medical Air	112-42	\$249.00
DISS Instrument Air	277-42	\$349.00
DISS Nitrogen	140-42	\$349.00



# HIGH PRESSURE BULK HOSE SOLD BY THE FOOT



- \* Available in custom lengths and configurations
- \* Available in a variety of gasses
- \* Sold by The Foot \* Built to your needs\*
- \* 200 PSI Max

GAS	CONDUCTIVE	NON-CONDUCTIVE
OXYGEN	22-01	20-01
MEDICAL AIR	23-01	21-01
VACUUM 1/4"	25-01	
VACUUM 5/16"	32-01	
N2O	24-01	
CO2	28-01	
WAGD 1/4"	24-01	
WAGD 5/16"	33-01	
NITROGEN	44-01	
INSTRUMENT AIR	34-01	
GENERAL PURPOSE	35-01	

**\*\*Hose Lengths Sold By The Foot\*\***

**FERRULE FOR 1/4" HOSE**



Y624A      10/PACK

**HEAVY DUTY HOSE RETRACTOR**



60-64

**HEAVY DUTY SECURITY CABLE IDEAL FOR SUCTION REGULATORS**



#OH-167      \$5.50 Each

## STANDARD HIGH PRESSURE HOSE

10 FT NON CONDUCTIVE RESPIRATORY \$54.00		
#100-10-20-100	Green	Oxygen
#110-10-21-110	Yellow	Air
10 FT CONDUCTIVE OR/ANESTHESIA \$89.00		
#100-10-22-100	Green	Oxygen
#110-10-23-110	Yellow	Air
#120-10-24-120	Blue	Nitrous Oxide
#130-10-25-130	White	Vacuum
140-10-26-140	Black	Nitrogen
#150-10-27-150	Lavender	EVAC
#180-10-28-180	Gray	Carbon Dioxide

**DISS Female x DISS Female, 10 FT Long**





**13" Chord/Hose Wraps \$59.00/10 PACK**




Red	WTV13-R
Blue	WTV13-BL
Black	WTV-13-BK

## CHEMETRON FEMALE COUPLERS

1/4" HOSE BARB		
		
CHEMETRON / NCG		
400-77	OXYGEN	\$69.00
410-77	AIR	\$69.00
430-77	VACUUM	\$79.00
420-77	N2O	\$89.00
450-77	WAGD	\$149.00
480-77	CO2	\$149.00


DISS MALE		
		
CHEMETRON / NCG		
400-55	OXYGEN	\$69.00
410-55	AIR	\$69.00
430-55	VACUUM	\$79.00
420-55	N2O	\$89.00
450-55	WAGD	\$149.00
480-55	CO2	\$149.00


DISS FEMALE		
		
CHEMETRON / NCG		
55-400	OXYGEN	\$69.00
55-410	AIR	\$69.00
55-430	VACUUM	\$79.00
55-420	N2O	\$89.00
55-450	WAGD	\$139.00
55-480	CO2	\$149.00

1/4" NPT MALE		
		
CHEMETRON / NCG		
400-14	OXYGEN	\$65.00
410-14	AIR	\$65.00
430-14	VACUUM	\$79.00
420-14	N2O	\$89.00
450-14	WAGD	\$139.00
480-14	CO2	\$139.00

## PURITAN BENNETT FEMALE COUPLERS

1/4" HOSE BARB		
		
Puritan Bennett		
200-77	OXYGEN	\$69.00
210-77	AIR	\$69.00
230-77	VACUUM	\$79.00
220-77	N2O	\$89.00
250-77	WAGD	\$139.00
280-77	CO2	\$139.00

DISS MALE		
		
Puritan Bennett		
200-55	OXYGEN	\$74.00
210-55	AIR	\$74.00
230-55	VACUUM	\$84.00
220-55	N2O	\$96.00
250-55	WAGD	\$139.00
280-55	CO2	\$139.00

DISS FEMALE		
		
Puritan Bennett		
55-200	OXYGEN	\$74.00
55-210	AIR	\$74.00
55-230	VACUUM	\$84.00
55-220	N2O	\$96.00
55-250	WAGD	\$139.00
55-280	CO2	\$139.00

1/4" NPT MALE		
		
Puritan Bennett		
200-14	OXYGEN	\$69.00
210-14	AIR	\$69.00
230-14	VACUUM	\$79.00
220-14	N2O	\$89.00
250-14	WAGD	\$139.00
280-14	CO2	\$139.00

## Custom High Pressure Hoses



### 10" ADJUSTABLE WRENCH FOR TANK REGULATORS

10-AW	\$28.00
-------	---------

### 1/4" HOSE BARB CHEMETRON



PART #	GAS	PRICE
401-77	O2	\$29.00
411-77	AIR	\$29.00
421-77	N2O	\$46.00
431-77	VAC	\$36.00
451-77	EVAC	\$89.00
481-77	CO2	\$89.00

### 1/4" HOSE BARB PB



PART #	GAS	PRICE
201-77	O2	\$27.00
211-77	AIR	\$32.50
221-77	N2O	\$44.00
231-77	VAC	\$36.00
251-77	EVAC	\$69.00
281-77	CO2	\$89.00

## MEDSTAR FEMALE COUPLERS

1/4" HOSE BARB		
MEDSTAR		
300-77	OXYGEN	\$79.00
310-77	AIR	\$79.00
330-77	VACUUM	\$89.00
320-77	N2O	\$89.00
353-77	WAGD	\$149.00
	CO2	\$159.00

DISS MALE		
MEDSTAR		
300-55	OXYGEN	\$79.00
310-55	AIR	\$79.00
330-55	VACUUM	\$79.00
320-55	N2O	\$139.00
352-55	WAGD	\$149.00

DISS FEMALE		
MEDSTAR		
55-300	OXYGEN	\$89.00
55-310	AIR	\$89.00
55-330	VACUUM	\$89.00
55-320	N2O	\$149.00
55-352	WAGD	\$159.00

1/4" NPT MALE		
MEDSTAR		
300-14	OXYGEN	\$79.00
310-14	AIR	\$79.00
330-14	VACUUM	\$89.00
320-14	N2O	\$89.00
352-14	WAGD	\$149.00
	CO2	\$159.00

## OHIO FEMALE COUPLERS


1/4" HOSE BARB		
OHIO		
500-77	OXYGEN	\$79.00
510-77	AIR	\$89.00
530-77	VACUUM	\$89.00
520-77	N2O	\$96.00
550-77	WAGD	\$139.00
580-77	CO2	\$139.00


DISS MALE		
OHIO		
500-55	OXYGEN	\$79.00
510-55	AIR	\$89.00
530-55	VACUUM	\$89.00
520-55	N2O	\$96.00
550-55	WAGD	\$139.00
580-55	CO2	\$139.00

DISS FEMALE		
OHIO		
55-500	OXYGEN	\$79.00
55-510	AIR	\$89.00
55-530	VACUUM	\$96.00
55-520	N2O	\$96.00
55-550	WAGD	\$139.00
55-580	CO2	\$139.00

1/4" NPT MALE		
OHIO		
500-14	OXYGEN	\$75.00
510-14	AIR	\$85.00
530-14	VACUUM	\$89.00
520-14	N2O	\$96.00
550-14	WAGD	\$139.00
580-14	CO2	\$139.00

## Custom High Pressure Hoses

1/4" HOSE BARB MEDSTAR		
		
PART #	GAS	PRICE
302-77	O2	\$54.50
312-77	AIR	\$54.50
322-77	N2O	\$69.00
332-77	VAC	\$59.00
354-77	EVAC	\$94.00
355-77	CO2	\$94.00


1/4" HOSE BARB OHIO		
		
PART #	GAS	PRICE
501-77	O2	\$29.50
511-77	AIR	\$29.50
521-77	N2O	\$46.00
531-77	VAC	\$36.50
551-77	EVAC	\$69.00
581-77	CO2	\$89.00




**\*\* ASK ABOUT OUR \*\*  
90° SWIVEL HOSES**

**DISS BODY ADAPTERS**

DISS MALE BY NPT MALE		
1/8" PART#	1/4" PART#	GAS
102-18	102-14	OXYGEN
112-18	112-14	AIR
122-18	122-14	NITROUS OXIDE
132-18	132-14	VACUUM
142-18	142-14	NITROGEN
Not Available	152-14	EVAC
Not Available	182-14	CARBON DIOXIDE
Not Available	197-14	HELIUM
162-18	162-14	INSTRUMENT AIR




1/8" NPT MALE




1/4" NPT MALE

DISS MALE BY NPT MALE WITH CHECK VALVE		
1/8" PART#	1/4" PART#	GAS
102-CV-18	102-CV-14	OXYGEN
112-CV-18	112-CV-14	AIR
122-CV-18	122-CV-14	NITROUS OXIDE
Not Available	132-CV-14	VACUUM
142-CV-18	142-CV-14	NITROGEN
Not Available	152-CV-14	CARBON DIOXIDE




1/4" NPT MALE




1/4" NPT MALE

DISS MALE BY NPT FEMALE		
1/8" PART#	1/4" PART#	GAS
102-18F	102-14F	OXYGEN
112-18F	112-14F	AIR
122-18F	122-14F	NITROUS OXIDE
132-18F	132-14F	VACUUM
142-18F	142-14F	NITROGEN




1/4" NPT FEMALE

DISS MALE BY NPT FEMALE	
PART#	GAS
102-77-14	OXYGEN
110-77-14	AIR
122-77-14	NITROUS OXIDE
132-77-14	VACUUM
142-77-14	NITROGEN



1/4" HOSE BARB

DISS MALE BY DISS MALE	
PART#	DESCRIPTION
102-102	Oxygen - DISS 1240, 9/16"




DISS MALE




**DISS BODY ADAPTERS**

<b>DISS FEMALE HEX NUT BY NPT</b>		
<b>1/8" PART#</b>	<b>1/4" PART#</b>	<b>GAS</b>
101-18	101-14	OXYGEN
111-18	111-14	AIR
121-18	121-14	NITROUS OXIDE
131-18	131-14	VACUUM
141-18	141-14	NITROGEN
191-18	181-14	CARBON DIOXIDE
161-18	161-14	INSTRUMENT AIR
151-18	151-14	EVAC




**HEX NUT BY NPT**

<b>DISS FEMALE HAND TIGHT BY NPT</b>		
<b>1/8" PART#</b>	<b>1/4" PART#</b>	<b>GAS</b>
100-18	100-14	OXYGEN
110-18	110-14	AIR
120-18	120-14	NITROUS OXIDE
130-18	130-14	VACUUM
140-18	140-14	NITROGEN
180-18	180-14	CARBON DIOXIDE
150-18	150-14	EVAC
160-18	160-14	INSTRUMENT AIR




**HAND TIGHT BY NPT**

<b>DISS FEMALE HEX NUT BY 1/4" HOSE BARB</b>	
<b>PART#</b>	<b>GAS</b>
101-77-14	OXYGEN
111-77-14	AIR
121-77-14	NITROUS OXIDE
131-77-14	VACUUM
141-77-14	NITROGEN
151-77-14	EVAC
181-77-14	CARBON DIOXIDE
198-77-14	HELIUM
161-77-14	INSTRUMENT AIR







**HEX NUT BY 1/4" HOSE BARB**


<b>DISS FEMALE HAND TIGHT BY 1/4" HOSE BARB</b>	
<b>PART#</b>	<b>GAS</b>
100-77-14	OXYGEN
110-77-14	AIR
120-77-14	NITROUS OXIDE
130-77-14	VACUUM
140-77-14	NITROGEN
150-77-14	EVAC
180-77-14	CARBON DIOXIDE
190-77-14	HELIUM
160-77-14	INSTRUMENT AIR





**HAND TIGHT BY 1/4" HOSE BARB**


**NPT FITTINGS 200 PSI MAX / TUBING NIPPLES**


HEX NIPPLE		CONNECTOR		COUPLER		90° STREET ELBOW	
							
PART#	DESCRIPTION	PART#	DESCRIPTION	PART#	DESCRIPTION	PART#	DESCRIPTION
MH18-18	1/8" M NPT x 1/8" M NPT	FC18-M18	1/8" F NPT x 1/8" M NPT	F18-F18	1/8" F NPT x 1/8" F NPT	FE18-M18	1/8" F NPT x 1/8" M NPT
MH14-14	1/4" M NPT x 1/4" M NPT	FC18-M14	1/8" F NPT x 1/4" M NPT	F14-F14	1/4" F NPT x 1/4" F NPT	FE14-M18	1/4" F NPT x 1/8" M NPT
MH18-14	1/8" M NPT x 1/4" M NPT	FC14-M18	1/4" F NPT x 1/8" M NPT	F14-F18	1/4" F NPT x 1/8" F NPT	FE14-M14	1/4" F NPT x 1/4" M NPT
		FC14-M14	1/4" F NPT x 1/4" M NPT				


90° ELBOW	
	
PART#	DESCRIPTION
18-ENPT	1/8" F NPT
14-ENPT	1/4" F NPT


WYE BLOCK	
	
PART#	DESCRIPTION
22-18-18	1/8" All Ports
22-18-14	1/8" Inlet, 1/4" Outlets


FEMALE TEE	
	
PART#	DESCRIPTION
T-18-18	1/8" F NPT
T-14-14	1/4" F NPT


POWER TAKE OFF	
	
PART#	DESCRIPTION
P-24	OXYGEN
P-16	MEDICAL AIR

1/4" ID FEMALE NPT TO HOSE BARB	
	
PART#	DESCRIPTION
F-18-77	1/8" NPT Female
F-14-77	1/4" NPT Female

1/4" ID MALE NPT TO HOSE BARB	
	
PART#	DESCRIPTION
M-18-77	1/8" NPT Male
M-14-77	1/4" NPT Male

SHUT OFF VALVE	
	
PR-7216	\$79.00

90 Degree Swivel Adaptor to Hose Barb	
	
PART#	DESCRIPTION
FPT-2490S	1/8" NPT Female X 1/4" I.D. Hose Barb

DISS Female X 1/8" NPT Male- 1 1/2" Nipple Length	
	
PART#	DESCRIPTION
1249-2	1/8" NPT Female

YOKE TEE

013-00
\$7.50 Each

Wing Nut Nipple 0158

Nipple W/Chain 0158C




SURE SEAL WASHERS		
7-53	Aluminum, Viton	Pack of 50

**MEDICAL GRADE BALL VALVES/  
LOKRINGS/CHECK VALVES**



**BALL VALVES**



BALL VALVES		DIMENSIONS (Inches)						
Valve Size	Part Number	A	B	D	E	F	G	H
1/2"	41050-01	2.62	3.35	1.69	7.87	6.37	21.00	5.06
3/4"	41075-01	2.93	4.33	1.89	8.25	6.06	21.00	4.60
1"	41100-01	3.44	4.33	2.05	8.58	5.87	21.00	4.15
1-1/4"	41125-01	3.96	5.12	2.56	8.75	5.50	21.00	3.52
1-1/2"	41150-01	4.43	5.12	2.72	8.87	5.37	21.00	3.16
2"	41200-01	5.55	6.30	3.15	10.87	5.00	21.00	2.23
2-1/2"	41250-01	6.73	9.06	4.29	16.00	4.25	32.00	0.89
3"	41300-01	7.56	9.06	4.61	16.00	4.50	32.00	0.72
4"	41400-01	9.80	11.73	7.90	16.00	7.90	32.00	3.00

CHECK VALVE SIZES	
Check-Valve Size	Part Number
3/4"	51075-01
1"	51100-01
1-1/4"	51125-01
1-1/2"	51150-01
2"	51200-01
2-1/2"	51250-01
3"	51300-01
4"	51400-01

**NITROGEN CONTROL PANELS**

**AMICO**  
COMPLETE ASSEMBLY



DISS Male      N-CONP-E-REL  
Schrader Female    N-CONS-E-REL


**Nitrogen Control Retrofit Ohmeda**



Replacement Front Panel Only  
DISS Male      N-CONP-E-OHM  
Schrader Female    N-CONS-E-OHM


**BEACON**

**Nitrogen Control Panel II**




Complete Assembly  
DISS Male      6-120267-01  
Schrader Female    6-120277-0S

**Nitrogen Control Retrofit Ohmeda**




Replacement Front Panel Only  
DISS Male      6-120901-00  
Schrader Female    6-120901-0S

**Nitrogen Control Retrofit Puritan Bennett**




Replacement Front Panel Only  
DISS Male      6-120900-00  
Schrader Female    6-120900-0S

**Nitrogen Control for Ceiling Column**



Replacement Front Panel Only  
DISS Male      6-120396-00  
Schrader Female    6-120396-0S

**ALLIED**



COMPLETE ASSEMBLY

DISS Male      75-20-0005  
Schrader Female    75-20-0007

Replacement Front Panel Only

DISS Male      75-20-0015  
Schrader Female    75-20-0017

# GAS ANALYSIS & OUTLET VERIFICATION

**\*\*Member of the NFPA\*\***

**NFPA 99 requires that after any renovation, modification or addition, a medical gas or vacuum system is to be inspected in accordance with the standards set by the NFPA.**

**\*\*NFPA 99 & NFPA 55 REQUIRE ANNUAL OR SEMI-ANNUAL\*\*  
**\*\*INSPECTIONS OF PATIENT CARE FACILITIES\*\*****

## ValueMed Offers Many Different Testing Options

\*Medical Gas Outlets- Function & Operation

- Oxygen
- Medical Air
- Nitrous Oxide
- Nitrogen
- Carbon Dioxide
- Instrument Air
- Vacuum
- WAGD



\*Zone Valves-Function & Operation

\*Alarm Panels-Area and Master-Function & Operation

\*Medical Gas Manifolds-Function & Operation

\*Medical Air Compressor-Function & Operation

\*Medical Vacuum Pump-Function & Operation

\*Medical Gas Purity Analysis

\*Dew Point Determination

\*Trace Gas Analysis


\*Lab Hoods-Function & Operation

\*Carbon Monoxide Safety Tests

\*Methane Safety Tests

\*Trace Hydrocarbons Safety Tests

\*Employee Exposure Monitoring



Current Code Book  
PART # NFPA-HCFC



## Full line of Shutdown Equipment

- Carts-Regulators-Suction Pumps-
- Oxygen-Air-Nitrous Oxide-

**STATE OF THE ART EQUIPMENT  
& QUALIFIED TECHNICIANS**



MBIV

\*Reports Prepared Instantly\*

\*Purity Results On Site\*





**BEACON MEDAES SERIES "B" WALL OUTLETS**



Item	Concealed Wall Outlet				
	Gas Service	Color Code	Complete Assembly	Rough-in Assembly	Latch-Valve Assembly
Geometric Quick-Connect (U.S.)	O <sub>2</sub>	White	<input type="checkbox"/> 6-121020-00	<input type="checkbox"/> 6-233110-00	<input type="checkbox"/> 6-230930-00
	N <sub>2</sub> O	Blue	<input type="checkbox"/> 6-121021-00	<input type="checkbox"/> 6-233111-00	<input type="checkbox"/> 6-230931-00
	AIR	Yellow	<input type="checkbox"/> 6-121022-00	<input type="checkbox"/> 6-233112-00	<input type="checkbox"/> 6-230932-00
	VAC	White	<input type="checkbox"/> 6-121023-00	<input type="checkbox"/> 6-233113-00	<input type="checkbox"/> 6-230933-00
	WAGD	Purple	<input type="checkbox"/> 6-121029-00	<input type="checkbox"/> 6-233119-00	<input type="checkbox"/> 6-230939-00

Add "P" to part number for PURITIAN BENNETT  
 Add "C" to part number for CHEMETRON  
 Add "O" to part number for OHIO  
 Add "M" to part number for MEDSTAR/OXEQUIP

**BEACON DISS OUTLETS**

Item	Concealed Wall Outlet				
	Gas Service	Color Code	Complete Assembly	Rough-in Assembly	Latch-Valve Assembly
Series B DISS (U.S.)	O <sub>2</sub>	White	<input type="checkbox"/> 6-121100-00	<input type="checkbox"/> 6-233110-00	<input type="checkbox"/> 6-230910-00
	N <sub>2</sub> O	Blue	<input type="checkbox"/> 6-121101-00	<input type="checkbox"/> 6-233111-00	<input type="checkbox"/> 6-230911-00
	AIR	Yellow	<input type="checkbox"/> 6-121102-00	<input type="checkbox"/> 6-233112-00	<input type="checkbox"/> 6-230912-00
	VAC	White	<input type="checkbox"/> 6-121103-00	<input type="checkbox"/> 6-233113-00	<input type="checkbox"/> 6-230913-00
	N <sub>2</sub>	Black	<input type="checkbox"/> 6-121104-00	<input type="checkbox"/> 6-233114-00	<input type="checkbox"/> 6-230914-00
	Instrument Air	Red	<input type="checkbox"/> 6-121108-00	<input type="checkbox"/> 6-233118-00	<input type="checkbox"/> 6-230916-00
	WAGD	Purple	<input type="checkbox"/> 6-121109-00	<input type="checkbox"/> 6-233119-00	<input type="checkbox"/> 6-230919-00
	CO <sub>2</sub>	Gray	<input type="checkbox"/> 6-121110-00	<input type="checkbox"/> 6-233120-00	<input type="checkbox"/> 6-230920-00

Standard Trim Plate 5" #6-325161-00	Console Trim Plate 3.25" #6-325160-00
-------------------------------------	---------------------------------------

**A COMPLETE OUTLET ASSEMBLY NEEDS 3 PART NUMBERS  
 (ROUGH-IN ASSEMBLY, LATCH-VALVE ASSEMBLY & TRIM PLATE)**

**OHIO WALL OUTLETS**


**OHIO DIAMOND III / GEMINI**



Concealed Wall Outlet				
Item	Gas Service	Color Code	Rough-In Assembly	Finish Assembly
<b>Wall/Ceiling Quick Connect</b>	O <sub>2</sub>	Green	<input type="checkbox"/> 234510-00	<input type="checkbox"/> 234610-00
	N <sub>2</sub> O	Blue	<input type="checkbox"/> 234511-00	<input type="checkbox"/> 234611-00
	AIR	Yellow	<input type="checkbox"/> 234512-00	<input type="checkbox"/> 234612-00
	VAC	White	<input type="checkbox"/> 234513-00	<input type="checkbox"/> 234613-00
	WAGD	Purple	<input type="checkbox"/> 234519-00	<input type="checkbox"/> 234619-00
<b>Console Quick-Connect</b>	O <sub>2</sub>	Green	<input type="checkbox"/> 234410-00	<input type="checkbox"/> 234810-00
	N <sub>2</sub> O	Blue	<input type="checkbox"/> 234411-00	<input type="checkbox"/> 234811-00
	AIR	Yellow	<input type="checkbox"/> 234412-00	<input type="checkbox"/> 234812-00
	VAC	White	<input type="checkbox"/> 234413-00	<input type="checkbox"/> 234813-00
	WAGD	Purple	<input type="checkbox"/> 234419-00	<input type="checkbox"/> 234819-00
<b>Wall/Ceiling DISS</b>	O <sub>2</sub>	Green	<input type="checkbox"/> 234510-00	<input type="checkbox"/> 234710-00
	N <sub>2</sub> O	Blue	<input type="checkbox"/> 234511-00	<input type="checkbox"/> 234711-00
	AIR	Yellow	<input type="checkbox"/> 234512-00	<input type="checkbox"/> 234712-00
	VAC	White	<input type="checkbox"/> 234513-00	<input type="checkbox"/> 234713-00
	N <sub>2</sub>	Black	<input type="checkbox"/> 234514-00	<input type="checkbox"/> 234714-00
	WAGD	Purple	<input type="checkbox"/> 234519-00	<input type="checkbox"/> 234719-00
	CO <sub>2</sub>	Gray	<input type="checkbox"/> 234520-00	<input type="checkbox"/> 234720-00
<b>Accessories</b>	Vacuum Slide		<input type="checkbox"/> 234500-00	<input type="checkbox"/> 234600-00
	Blanket Outlet		<input type="checkbox"/> 234500-00	<input type="checkbox"/> 234601-00

Medaes\* and Diamond\* are trademarks of Hillenbrand Industries  
Ohio\* and Ohmeda\* are trademarks of Datex-Ohmeda

**OHIO DIAMOND CARE**

	FINISH ASSEMBLY ONLY		ROUGH-IN ASSEMBLY ONLY	
	OXYGEN	6803-8121-800	OXYGEN	6803-8139-800
	MEDICAL AIR	6803-8121-803	MEDICAL AIR	6803-8139-803
	NITROUS OXIDE	6803-8121-801	NITROUS OXIDE	6803-8139-801
	VACUUM	6803-8121-802	VACUUM	6803-8139-802
	Add "C" to part number for CHEMETRON			
Add "P" to part number for PURITAN BENNETT				
Add "D" to part number for DISS				
Add "CL" to part number for CONSOLE MOUNT				

**GAS ANALYSIS & OUTLET VERIFICATION SERVICES**  
**MEDICAL GAS & VACUUM OUTLETS - ALARMS - VALVE BOXES**



- \* Over 30 Years of Testing Experience \*
- \* Complete Documentation for Permanent Records \*
- \* Annual Testing \*
- \* All Testing Performed in Accordance with NFPA 99C \*






**Value Med**


**TO PLACE AN ORDER**


**FAX: 888-676-0711**


**EMAIL: [valuemed.sales@gmail.com](mailto:valuemed.sales@gmail.com)**

### CHEMETRON / MEDSTAR WALL OUTLETS

	CHEMETRON 500 SERIES			
	COMPLETE OUTLET ASSEMBLY		VALVE BODY ASSEMBLY	COVER PLATE
	OXYGEN	64-01-5001	64-90-0620	64-91-1501
	VACUUM	64-01-5002	64-90-0621	64-91-1502
	MEDICAL AIR	64-01-5003	64-90-0620	64-91-1503
	NITROUS OXIDE	64-01-5004	64-90-0620	64-91-1504
	WAGD	64-01-5006	64-90-0621	64-91-1506
	CARBON DIOXIDE	64-01-5007	64-90-0620	64-91-1507
	O2/CO2	64-01-5008	64-90-0620	64-91-1508
	SLIDE BRACKET	64-06-001		




	DISS OUTLETS 500 SERIES		
	COMPLETE OUTLET ASSEMBLY	REPLACEMENT VALVE KIT	
	OXYGEN	64-04-5601	64-90-0631
	MEDICAL AIR	64-04-5603	64-90-0634
	NITROUS OXIDE	64-04-5604	64-90-0633
	VACUUM	64-04-5602	64-90-0632
	WAGD	64-04-5606	64-90-0638
	N2	64-04-5607	64-90-0636
	CO2	64-04-5608	64-90-0635

	CHEMETRON 400 SERIES			
	COMPLETE OUTLET ASSEMBLY	REPLACEMENT VALVE KIT	COVER PLATE	
	OXYGEN	64-20-001	64-90-0180	64-92-0176
	MEDICAL AIR	64-20-002	64-90-0180	64-92-0178
	NITROUS OXIDE	64-20-006	64-90-0180	64-91-4003
VACUUM	64-20-005	64-90-0181	64-92-0177	

	MEDSTAR		
	COMPLETE OUTLET	REPLACEMENT VALVE KIT	
	OXYGEN	1930-2111	1227-X
	MEDICAL AIR	1930-2411	1227-X
	NITROUS OXIDE	1930-2311	1227-X
	VACUUM	1930-2211	1227-X
	WAGD	1930-2811	1227-X
	CARBON DIOXIDE	1930-2511	1227-X
SLIDE BRACKET	1857-A		



**PRESSURE AND VACUUM GAUGES \$39.00**

		Gauge with Rubber Protector 
MG-C2	MG-C23	MJ-C1P





MG-C2	4000 PSI	2" OD	1/4" NPT	Chrome
MG-C3	100 PSI	1 1/2" OD	1/8" NPT	Chrome
MG-C4	100 PSI	2" OD	1/4" NPT	Chrome
MG-C7	15 LPM	2" OD	1/4" NPT	Chrome
MG-C23	4000 PSI	1 1/2" OD	1/8" NPT	Chrome
MG-C75	3000 PSI	1 1/2" OD	1 1/2" NPT	Chrome
MG-C75P	3000 PSI	1 1/2" OD	1 1/2" NPT	Chrome

	MG-C29		MG-C58		MG-C14
			Gomco 300 Gauge 01-90-3059		

**Leak Test Solution**



LD-8





BOTTOM PORT					BACK PORT				
<b>MG-C9</b>	4000 psi	2" OD	1/4" NPT	Chrome	<b>MG-B57</b>	100 psi	1-1/2" OD	1/8" NPT	Black
<b>MG-C11</b>	100 psi	2" OD	1/4" NPT	Chrome	<b>MG-B59</b>	300 psi	1-1/2" OD	1/8" NPT	Black
<b>MG-C13</b>	4000 psi	1-1/2" OD	1/8" NPT	Chrome	<b>Replacement Lenses</b>				
<b>MG-C14</b>	30 psi	2" OD	1/4" NPT	Chrome	<b>MGL-15</b>	Screw-On Type for 1-1/2" Chrome Gauges			
<b>MG-C15</b>	60 psi	2" OD	1/4" NPT	Chrome	<b>MGL-15WS</b>	Snap-In Type for 1-1/2" Gauges			
<b>MG-C16</b>	200 psi	2" OD	1/4" NPT	Chrome	<b>MGL-20</b>	Screw-On Type for 2" Chrome and Black Gauges			
<b>MG-C17</b>	400 psi	2" OD	1/4" NPT	Chrome	<b>MGL-70</b>	Snap-In Type for 2-1/2" Suction Gauges			
<b>MG-C18</b>	600 psi	2" OD	1/4" NPT	Chrome					
<b>MG-C19</b>	1000 psi	2" OD	1/4" NPT	Chrome					
NITROUS OXIDE Blue face w/white letters					VACUUM				
<b>MG-C29</b>	2000 psi	2" OD	1/4" NPT	Chrome	<b>MG-C58</b>	30 inHg Vac	1-1/2" OD	1/8" NPT	Chrome
<b>MG-C30</b>	100 psi	2" OD	1/4" NPT	Chrome	<b>MG-B76</b>	200 mm HG/	2" OD	1/8" NPT	Black
					Bottom Port				

**Teflon Tape**



TFE-40

**REPLACEMENT SUCTION REGULATOR GAUGES \$89.00**

			
BOEHRINGER 1/8" NPT	OHIO/ OHMEDA	PRECISION	CHEMETRON
BG-300	HG-300	PG-300	CG-300
OXEQUIP 1/4" NPT	AMVEX		
OX-300	AM-300		



**OUTLET ACCESSORIES / WALL OUTLET TESTER**



**Hospital Wall Outlet Test Gauge for Checking Gas Line Flow & Pressure**  
ALL \$369.00

OUTLET	O2	AIR	N2O	VAC	N2	INST AIR
DISS	20-DT	21-DT	24-DT	25-DT	26-DT	34-DT
OHIO	20-OT	21-OT	24-OT	25-OT	26-OT	34-OT
CHEMETRON	20-CT	21-CT	24-CT	25-CT	26-CT	34-CT
MEDSTAR	20-MT	21-MT	24-MT	25-MT	26-MT	34-MT
PURITAN BENNETT	20-PT	21-PT	24-PT	25-PT	26-PT	34-PT



MASTER ALARM PRESSURE / VACUUM SWITCH		
DEMAND CHECK FITTINGS (GAS SPECIFIC)		
	ALLIED	BEACON
OXYGEN	13-0100	6-129030-10
MEDICAL AIR	13-0103	6-129030-12
NITROUS OXIDE	13-0102	6-129030-11
CARBON DIOXIDE	13-0104	6-129030-20
VACUUM	13-0107	6-129030-13
WAGD	13-0107	6-129030-19
NITROGEN	*13-0105	6-129030-14
INSTRUMENT AIR	*13-0130	6-129030-18
*Requires 2 Switches		

**O-Ring Lubricant**



64-90-2111

For Medical Gas & Vacuum Outlets


**ENVIROVAC**



64-90-VAC

Outlet Cleaner

**Zone Valve Box Replacement Gauges**

1.5" Gauge		
1/8" NPT Back Port		
100 PSI	MG-BW-70	
300 PSI	MG-BW-71	
30 In/Hg	MG-BW-72	

**4 1/2" Gauge With 1/4" NPT**



Steel-Case Gauges	
Part Number	Item
<input type="checkbox"/> 6-130111-00	Vacuum Gauge (30 in. Hg)
<input type="checkbox"/> 6-130112-00	Low Pressure Gauge (0-100 psi)*
<input type="checkbox"/> 6-130113-00	High Pressure Gauge (0-300 psi)*



Sold as a 2 Piece Kit



**MASTER ALARM**

1/4" or 1/8" NPT x DISS Gas Fitting with Check Valve (GAS SPECIFIC)

OXYGEN	MA-102-101
MEDICAL AIR	MA-112-111
NITROUS OXIDE	MA-122-121
NITROGEN	MA-142-142
CARBON DIOXIDE	MA-182-181
INSTRUMENT AIR	MA-162-161
VACUUM	MA-132-131
WAGD	MA-152-151



**YOUR MEDICAL HARDWARE SPECIALISTS**

**24 HOUR SERVICE**

**888-998-VMED (8633)**

**EQUIPMENT REPAIR SERVICES**

**\*\*CALL US TODAY FOR EQUIPMENT REPAIR PICK-UP\*\***

Medical Gas  
Regulators & Flowmeters



Oxygen  
Medical Air  
Heli-Ox



Carbon Dioxide  
Nitrogen  
Nitrous Oxide

Suction Regulators



Repair or Upgrade Your  
Overflow Safety Trap

**\$39.50**

**Call For  
Estimates**



\* Chemetron \* Ohmeda \* \* Precision Medical \* Boehringer \*

Air/Oxygen Blenders

Overhauls  
Water Traps  
Hoses  
Brackets



Bird - Bio Med Devices - Precision Medical  
*Recommended Overhall every 2 years*

**FULL LINE OF CUSTOM HOSES**

**AVAILABLE**

\*High Pressure Hoses \*Hose  
Assemblies\*Surgery\*

**\*Anesthesia\*Respiratory\***

**Thank You Kevin McDermott**

**(708)369-3675**

**KMCDMED@AOL.COM**



Oxygen Analyzers & Monitors  
Replacement Sensors



MiniOx \* Hudson \*Maxtec



Ventilator Support Arms



Horizon Gas Block Sales and Service



Laryngoscope  
Blades & Handles

\* Replacement\*  
Lamps



Cufflators



Manometers

Respirometers



Wright -Boehringer-Haloscale





Percussors



General Physiotherapy  
Meds Systems

## RESPIRATORY EQUIPMENT

**\* RENTALS \* REPAIRS \* SERVICE \* SHUTDOWNS \***

-  **VENTILATORS**
-  **VISIONS**
-  **V 60's**
-  **VAPOTHERMS**



## SHUTDOWN SERVICES

**BLENDER / FISHER  
PAYKEL HEATER  
ON ROLLING STAND**



*LET US HELP WITH YOUR EMERGENCY AND PLANNED SHUTDOWNS*

**COMPRESSORS AND  
SUCTION PUMPS**



**CYLINDER CARTS**



*DUAL & SINGLE  
CARTS AVAILABLE*

**REGULATORS**



### **SHUTDOWN EQUIPMENT**

- \* REGULATORS/ PIGTAILS**
- \* CYLINDER CARTS**
- \* HOSES**
- \* AIR COMPRESSORS**
- \* SUCTION PUMPS**
- STAND & TABLE TOP**

**YOUR MEDICAL HARDWARE SPECIALISTS**

**24 HOUR SERVICE**

**888-998-VMED (8633)**



**Value Med**



**TO PLACE AN ORDER:**

**PHONE: 888-998-8633**

**FAX: 888-676-0711**

or

**EMAIL: [Valuemed.sales@gmail.com](mailto:Valuemed.sales@gmail.com)**

We are currently providing service and sales to: Froedtert and the Medical College of Wisconsin, Franciscan Alliance, Ascension, Indiana University, US Department of Veteran Affairs, Presence, University Of Wisconsin, NorthShore Health Systems, Aurora Healthcare, CHS, Cadence, Trinity, Wheaton Franciscan, Rush, Kindred, Community Health Systems, Northwestern Memorial Healthcare, Tenet Healthcare, Prohealth Care, Centegra Health System, Methodist Hospitals, University of Illinois Hospital and Health Sciences System, Advocate and the Adventist Group.

We would like to be your *Valued* partner in healthcare!

Our biomedical services are supplied to Aramark, Advocate Clinical Engineering, Sodexo, Philips Healthcare, Crothall, GE, Jannx Medical Systems and Trimedx and Ventech Medical.

